



200-Place High Profile Storage Box, Natural



Product Highlights

- 200-place polypropylene storage box with fixed 10x20 dividers for tubes up to 12.5mm, eg. screw cap microtubes.
- Hinged lid with positioning lugs for stable stacking
- Moulded reference grid inside
- Temperature use -90°C to 121°C
- Autoclavable
- Dimensions (WxDxH): 275 x 155 x 64 mm
- Max tube height: 53 mm

Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk




General Data

Art. No.	E2200-1000
Pack Size	1 Pcs. (1 Pack × 1 Pcs.)
Dimensions (W × D × H)	275 mm × 155 mm × 64 mm
Dimensions (W × D)	275 mm × 155 mm
Max. tube length	53 mm
External width	275 mm
Autoclavable	Yes
Material	Polypropylene
Color	Natural
Places	200

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk




More informations about 200-Place High Profile Storage Box, Natural

200-place high-profile polypropylene storage box for taller tubes, with fixed, 10x20 dividers for tubes up to 12.5mm, eg. screw cap microtubes.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	200-Place High Profile Storage Box, Natural Color: Natural	1 Pcs. (1 Pack × 1 Pcs.)	E2200-1000
	Create PDF Data Sheet		
	200-Place High Profile Storage Box, Black Color: Black	1 Pcs. (1 Pack × 1 Pcs.)	E2200-1008
	Create PDF Data Sheet		

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk