



## 75mm High Cardboard Box only, No Dividers, Red



### Product Highlights

- Cardboard storage box, no dividers
- Standard 133 x 133 mm footprint
- For use down to -80 °C
- Moisture-resistant coating
- Non-toxic water-based glue

## Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

Art. No.	A9625-0004
Pack Size	6 Pcs. (1 Pack × 6 Pcs.)
Dimensions (W × D × H)	133 mm × 133 mm × 75 mm
Dimensions (W × D)	133 mm × 133 mm
External width	133 mm
Material	Coated Cardboard
Color	Various
Places	16

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about 75mm High Cardboard Box only, No Dividers, Red

A range of cardboard storage boxes with a standard 133 x 133 mm footprint, available in four heights. Moisture-resistant coating makes the boxes a durable yet economical storage solution. Available in five variants to help you colour code your samples. Available with or without dividers.