



## 300 µl graduated TipOne® Tip, Natural, Racks (non-sterile), Case



### Product Highlights

- Graduated at 100 µl, 200 µl and 300 µl
- Universal tip fit
- High quality virgin polypropylene tips, with extremely smooth surfaces for reduced sample retention
- Certified RNase, DNase, DNA and Pyrogen free
- Racks can be retained and reused with TipOne® refills
- Rack base made with ~25 % recycled plastic
- All components, excluding tips, are recyclable
- Buy by the case to save on packaging and transportation. Reuse the case box to send back your recycling.

## Product Applications

Pipetting liquids

Dispensing liquids

Mixing liquids

Filling plates and reaction vessels

Phase extraction and removal of supernatants

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

Art. No.	S1110-9800-C
Pack Size	7,680 Tips (80 Racks × 96 Tips)
Color Code Wafer	Green
Sterile	No
Autoclavable	Yes
Filter tip	Non-filter
Volume range	1.0 - 300 µl
Volume	300 µl
Tip Type	standard
Material	Polypropylene
Graduations	Yes
Color	Natural
Sustainability attributes	<p>Racks can be reused with TipOne® refills.</p> <p>Rack base made with ~25 % recycled plastic.</p> <p>All components, excluding tips, are recyclable.</p> <p>Buy by the case to save on packaging and transportation.</p> <p>Reuse the case box to send back your recycling.</p>

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



For use with

**Abimed:** BioPette-E, Discovery Comfort, LABMATE Optima

**Biohit:** eLine, ePet, mLine, Picus, Proline

**BRAND:** Transferpette, Transferpette Electronic, Transferpette Pro, Transferpette S

**Capp:** CappAero Comfort, CappAero Multi, CappAero Universal, Capptronic Multi, Capptronic Single, ecopipette, Singlechannel

**Eppendorf:** Reference, Reference 2, Research, Research plus, Research pro, Xplorer

**Gilson:** Pipetman Concept, Pipetman L, Pipetman Neo, Pipetman P, Pipetman Ultra

**Hamilton:** SoftGrip Pipette, SofTouch

**HTL:** Clinipet+, Discovery Comfort, Discovery Pro, Labmate Pro, Optipette

**Nichiryo:** 00-NP7, 50F, 50V, NLE, NPLO, NPX

**Rainin:** Classic, EDP1, EDP-3 Plus, Pipet-Lite, Pipet-Plus

**Sartorius:** eLine, ePet, mLine, Picus, Picus NxT, Proline, Proline Plus, Tacta

**Socorex:** Acura 815 manual, Acura 825, Acura 855, Acura electro 925, Acura electro 955, Calibra digital

**STARLAB:** ErgoOne, ErgoOne E, MicroOne, StarPet, StarPet-E

**Thermo Fisher Scientific Inc.:** Finn Focus, Finn Focus Micro, Finn Focus Univ., Finn timer Digital, Finn timer F1/F2, Finn timer Novus, Finn timer Novus Turquoise, Finn timer Novus Yellow

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about 300 µl graduated TipOne® Tip, Natural, Racks (non-sterile), Case

High quality, reliable pipette tips with a universal tip fit – TipOne® has it all! If you are looking for a sustainable pipette tip system with accuracy and precision you can trust, TipOne® is the one tip for you.

Join our journey to sustainable pipetting

TipOne® tip racks are built to last. When you have used your tips, the rack can be reused again and again with the cost-effective TipOne® refills (sterile and non-sterile options available). TipOne® racks have been tested and they can withstand at least 100 autoclave cycles (rack with wafer inside) without damage.













The base of the rack is made with ~25 % recycled plastic. All components, excluding tips, are recyclable. Sign up for our TipOne® recycling service!

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)











## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">300 µl graduated TipOne® Tip, Natural, Racks (non-sterile)</a> <b>Sterile:</b> No <b>Pack Size:</b> 960 Tips (10 Racks × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (10 Racks × 96 Tips)	S1110-9800
 <a href="#">300 µl graduated TipOne® Tip, Natural, Racks (non-sterile), Case</a> <b>Sterile:</b> No <b>Pack Size:</b> 7,680 Tips (80 Racks × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	7,680 Tips (80 Racks × 96 Tips)	S1110-9800-C
 <a href="#">300 µl graduated TipOne® Tip, Natural, Racks (sterile)</a> <b>Sterile:</b> Yes <b>Pack Size:</b> 960 Tips (10 Racks × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (10 Racks × 96 Tips)	S1110-9810
 <a href="#">300 µl graduated TipOne® Tip, Natural, Racks (sterile), Case</a> <b>Sterile:</b> Yes <b>Pack Size:</b> 7,680 Tips (80 Racks × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	7,680 Tips (80 Racks × 96 Tips)	S1110-9810-C
 <a href="#">300 µl graduated TipOne® Tip, Natural, Refills (non-sterile)</a> <b>Sterile:</b> No <b>Pack Size:</b> 960 Tips (10 Refills × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (10 Refills × 96 Tips)	S1110-9700
 <a href="#">300 µl graduated TipOne® Tip, Natural, Refills (non-sterile), Case</a> <b>Sterile:</b> No <b>Pack Size:</b> 9,600 Tips (100 Refills × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	9,600 Tips (100 Refills × 96 Tips)	S1110-9700-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">300 µl graduated TipOne® Tip, Refills (sterile)</a> <b>Sterile:</b> Yes <b>Pack Size:</b> 960 Tips (10 Refills × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (10 Refills × 96 Tips)	S1110-9710
 <a href="#">300 µl graduated TipOne® Tip, Natural, Refills (sterile), Case</a> <b>Sterile:</b> Yes <b>Pack Size:</b> 7,680 Tips (80 Refills × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	7,680 Tips (80 Refills × 96 Tips)	S1110-9710-C
 <a href="#">300 µl graduated TipOne® Tip, Natural, Bag (non-sterile)</a> <b>Sterile:</b> No <b>Pack Size:</b> 1,000 Tips (1 Bag × 1000 Tips)  <a href="#">Create PDF Data Sheet</a>	1,000 Tips (1 Bag × 1000 Tips)	S1110-9000
 <a href="#">300 µl graduated TipOne® Tip, Natural, Bags (non-sterile), Case</a> <b>Sterile:</b> No <b>Pack Size:</b> 25,000 Tips (25 Bags × 1000 Tips)  <a href="#">Create PDF Data Sheet</a>	25,000 Tips (25 Bags × 1000 Tips)	S1110-9000-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)