



## StarStore 81 Zylar Resin Storage Box



### Product Highlights

- 81-place zylar resin storage box with fixed 9x9 dividers for tubes up to 13mm
- Transparent lift-off lid with frosted writing panel on sides
- Vent holes in base for safety and to aid lid removal
- Temperature use: -196 °C to 121 °C
- Dimensions (WxDxH): 130 x 130 x 45 mm
- Can be used in the SS Freezer Racks or the Aluminum Freezer.

## Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	E2301-5309
Material	Zylar resin
Pack Size	5 Pcs. (1 Pack × 5 Pcs.)
Places	81
Autoclavable	Yes
External width	130 mm
Dimensions (W × D × H)	130 mm × 130 mm × 45 mm
Dimensions (W × D)	130 mm × 130 mm
Color	Natural

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about StarStore 81 Zylar Resin Storage Box

81-place zylar resin storage box with fixed, 9x9 dividers for tubes up to 13mm, eg. 1.0ml cryogenic vials. For use in ultra-low temperatures. Standard 130x130mm footprint for use in standard freezer racking.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)