



10 μ l graduated TipOne® Tip, Natural, Racks (non-sterile)



Product Highlights

- Graduated at 2.5 μ l and 10 μ l
- Universal tip fit
- High quality virgin polypropylene tips, with extremely smooth surfaces for reduced sample retention
- Certified RNase, DNase, DNA and Pyrogen free
- Racks can be retained and reused with TipOne® refills
- Rack base made with ~25% recycled plastic
- All components, excluding tips, are recyclable

Product Applications

Pipetting liquids

Dispensing liquids

Mixing liquids

Filling plates and reaction vessels

Phase extraction and removal of supernatants

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de




General Data

Art. No.	S1111-3800
Pack Size	960 Tips (10 Racks × 96 Tips)
Color Code Wafer	Red
Sterile	No
Autoclavable	Yes
Filter tip	Non-filter
Volume range	0.1 - 10 µl
Volume	10 µl
Tip Type	standard
Material	Polypropylene
Graduations	Yes
Color	Natural
Sustainability attributes	Racks can be reused with TipOne® refills. Rack base made with ~25 % recycled plastic. All components, excluding tips, are recyclable.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



For use with

Abimed: BioPette-E, Discovery Comfort, LABMATE Optima

Biohit: eLine, ePet, mLine, Picus, Proline

BRAND: Transferpette Electronic, Transferpette Pro, Transferpette S

Capp: CappAero Comfort, CappAero Multi, CappAero Universal, CappAero Universal short, ecopipette, Singlechannel

Eppendorf: Reference, Reference 2, Research, Research plus, Research pro, Xplorer

Gilson: Pipetman Concept, Pipetman L, Pipetman Neo, Pipetman P, Pipetman Ultra

Hamilton: SoftGrip Pipette, SofTouch

HTL: Discovery Comfort, Discovery Pro, Labmate Pro, Optipette

Nichiryo: 00-NP7, 50F, NLE, NPLO, NPX

Rainin: Classic, EDP1, EDP-3 Plus, Pipet-Lite, Pipet-Plus

Sartorius: eLine, ePet, mLine, Picus, Picus NxT, Proline, Proline Plus, Tacta

Socorex: Acura 815 manual, Acura 825, Acura 855, Acura electro 925, Acura electro 955, Calibra digital

STARLAB: ErgoOne, ErgoOne E, MicroOne, StarPet, StarPet-E

Thermo Fisher Scientific Inc.: Finn Focus, Finn Focus Micro, Finnpipette Digital, Finnpipette F1/F2, Finnpipette Novus, Finnpipette Novus Pink

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about 10 µl graduated TipOne® Tip, Natural, Racks (non-sterile)

High quality, reliable pipette tips with a universal tip fit – TipOne® has it all! If you are looking for a sustainable pipette tip system with accuracy and precision you can trust, TipOne® is the one tip for you.













Join our journey to sustainable pipetting

TipOne® tip racks are built to last. When you have used your tips, the rack can be reused again and again with the cost-effective TipOne® refills (sterile and non-sterile options available). TipOne® racks have been tested and they can withstand at least 100 autoclave cycles (rack with wafer inside) without damage.

The base of the rack is made with ~25 % recycled plastic. All components, excluding tips, are recyclable. Sign up for our TipOne® recycling service!














All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 10 µl graduated TipOne® Tip, Natural, Racks (non-sterile) Sterile: No Pack Size: 960 Tips (10 Racks × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Racks × 96 Tips)	S1111-3800
 10 µl graduated TipOne® Tip, Natural, Racks (non-sterile), Case Sterile: No Pack Size: 7,680 Tips (80 Racks × 96 Tips)  Create PDF Data Sheet	7,680 Tips (80 Racks × 96 Tips)	S1111-3800-C
 10 µl graduated TipOne® Tip, Natural, Racks (sterile) Sterile: Yes Pack Size: 960 Tips (10 Racks × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Racks × 96 Tips)	S1111-3810
 10 µl graduated TipOne® Tip, Natural, Racks (sterile), Case Sterile: Yes Pack Size: 7,680 Tips (80 Racks × 96 Tips)  Create PDF Data Sheet	7,680 Tips (80 Racks × 96 Tips)	S1111-3810-C
 10 µl graduated TipOne® Tip, Natural, Refills (non-sterile) Sterile: No Pack Size: 960 Tips (10 Refills × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Refills × 96 Tips)	S1111-3700
 10 µl graduated TipOne® Tip, Natural, Refills (non-sterile), Case Sterile: No Pack Size: 9,600 Tips (100 Refills × 96 Tips)  Create PDF Data Sheet	9,600 Tips (100 Refills × 96 Tips)	S1111-3700-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 10 µl graduated TipOne® Tip, Natural, Stack rack (non-sterile) Sterile: No Pack Size: 960 Tips (1 Stack rack × 960 Tips)  Create PDF Data Sheet	960 Tips (1 Stack rack × 960 Tips)	S1111-3200
 10 µl graduated TipOne® Tip, Natural, Stack racks (non-sterile), Case Sterile: No Pack Size: 9,600 Tips (10 Stack racks × 960 Tips)  Create PDF Data Sheet	9,600 Tips (10 Stack racks × 960 Tips)	S1111-3200-C
 10 µl graduated TipOne® Tip, Natural, Stack rack (sterile) Sterile: Yes Pack Size: 960 Tips (1 Stack rack × 960 Tips)  Create PDF Data Sheet	960 Tips (1 Stack rack × 960 Tips)	S1111-3210
 10 µl graduated TipOne® Tip, Natural, Stack racks (sterile), Case Sterile: Yes Pack Size: 9,600 Tips (10 Stack racks × 960 Tips)  Create PDF Data Sheet	9,600 Tips (10 Stack racks × 960 Tips)	S1111-3210-C
 10 µl graduated TipOne® Tip, Natural, Bag (non-sterile) Sterile: No Pack Size: 1,000 Tips (1 Bag × 1000 Tips)  Create PDF Data Sheet	1,000 Tips (1 Bag × 1000 Tips)	S1111-3000
 10 µl graduated TipOne® Tip, Natural, Bags (non-sterile), Case Sterile: No Pack Size: 50,000 Tips (50 Bags × 1000 Tips)	50,000 Tips (50 Bags × 1000 Tips)	S1111-3000-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de