



## Centrifuge Tube Storage / Transport Box for 15ml Tubes, Green



### Product Highlights

- Green polypropylene storage / transparent box for centrifuge tubes
- Holds 36 x 15ml tubes
- Transparent, lift-off lid
- Base can be removed for use in waterbath
- Temperature use -90°C to 121°C
- Autoclavable
- Dimensions (WxDxH): 145 x 145 x 127 mm

## Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	E9025-3615
Dimensions (W × D × H)	145 mm × 145 mm × 127 mm
Dimensions (W × D)	145 mm × 145 mm
Autoclavable	Yes
Pack Size	2 Pcs. (1 Box × 2 Pcs.)
Rack Type	Tube Rack
Material	Polypropylene
Color	Green
Max. centrifuge tubes or screwcap microtubes	15ml tubes
Places	36

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about Centrifuge Tube Storage / Transport Box for 15ml Tubes, Green

Polypropylene storage box that holds 36 x 15ml tubes.





Tubes not included.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">Centrifuge Tube Storage / Transport Box for 50ml Tubes, Blue</a> <b>Places: 16</b>  <a href="#">Create PDF Data Sheet</a>	2 Pcs. (1 Box × 2 Pcs.)	E9025-1650
 <a href="#">Centrifuge Tube Storage / Transport Box for 15ml Tubes, Green</a> <b>Places: 36</b>  <a href="#">Create PDF Data Sheet</a>	2 Pcs. (1 Box × 2 Pcs.)	E9025-3615

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)