



Totally Tubular™, Mixed



Product Highlights

- Four polypropylene racks in a support tray.
- Ind. rack dimensions (WxDxH): 93 x 93 x 48 mm.

Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	I2344-1000
Pack Size	1 Piece (1 Rack × 1 Piece)
Dimensions (W × D × H)	93 mm × 93 mm × 48 mm
Dimensions (W × D)	93 mm × 93 mm
External width	93 mm
Autoclavable	Yes
Rack Type	Multi-Tube Rack
Material	Polypropylene
Color	Mixed
Max. centrifuge tubes or screwcap microtubes	2 x 50 ml tubes (30 mm) 6 x 15 ml tubes (17 mm) 9 x 1.5 / 2.0 ml tubes (12 mm) 12 x 0.5 ml tubes (11 mm)
Places	29

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about Totally Tubular™, Mixed





Comprises four, 29-place polypropylene racks in a support tray. Tubes not included.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	Totally Tubular™, Mixed Color: Mixed	1 Piece (1 Rack × 1 Piece)	I2344-1000
	 Create PDF Data Sheet		
	Totally Tubular™, Neon Mixed Color: Neon mixed	1 Piece (1 Rack × 1 Piece)	I2344-1040
	 Create PDF Data Sheet		

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de