



## Benchtop Cooler / Transport Case



### Product Highlights

- Holds 48 x 1.5/2.0ml tubes
- Handy note panel inside the lid
- Maintains temperatures of:
  - -15 °C or lower, for 4 hours if frozen to -20 °C
  - -15 °C or lower, for 6 hours if frozen to -23 °C
- Dimensions (WxDxH): 215 x 200 x 95 mm
- Weight: 2.5 kg

## Product Applications

Sample storage

Sample transportation

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

|                        |                           |
|------------------------|---------------------------|
| Art. No.               | I6348-1515                |
| Pack Size              | 1 Piece (1 Box × 1 Piece) |
| Dimensions (W × D × H) | 215 mm × 200 mm × 95 mm   |
| Dimensions (W × D)     | 215 mm × 200 mm           |
| External width         | 215 mm                    |
| Weight                 | 2.5 kg                    |
| Autoclavable           | No                        |
| Material               | Polypropylene             |
| Color                  | Black                     |
| Places                 | 48                        |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about Benchtop Cooler / Transport Case

Small, polypropylene, case-style holder with ice-gel tube holder insert for 48 x 1.5/2.0 ml tubes. Protects valuable enzymes from temperature fluctuation due to frostfree cycles, power loss or open freezer doors.

**NOTE:** Not autoclavable

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)