



Mini Cooler



Product Highlights

- Holds 8 x 1.5/2.0 ml tubes
- Hinged, transparent lid with snap-catch
- Maintains -15 °C or below for up to 1 hour
- Dimensions (WxDxH): 120 x 85 x 90 mm
- Weight: 750 g

Product Applications

Sample storage
Sample transportation

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	I6308-1515
Dimensions (W × D × H)	120 mm × 85 mm × 90 mm
Dimensions (W × D)	120 mm × 85 mm
Weight	700 g
Autoclavable	No
Pack Size	1 Piece (1 Rack × 1 Piece)
Color	Silver
Places	8

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about Mini Cooler

Small rack with ice-gel tube holder insert for 8 x 1.5/2.0 ml tubes. Perfect for the transportation of a small number of samples.

NOTE: Must be frozen to -20 °C for a minimum of 24 hours before use. Not autoclavable.