



## Mini Cooler



### Product Highlights

- Holds 8 x 1.5/2.0 ml tubes
- Hinged, transparent lid with snap-catch
- Maintains -15 °C or below for up to 1 hour
- Dimensions (WxDxH): 120 x 85 x 90 mm
- Weight: 750 g

## Product Applications

Sample storage  
Sample transportation

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

|                        |                            |
|------------------------|----------------------------|
| Art. No.               | I6308-1515                 |
| Dimensions (W × D × H) | 120 mm × 85 mm × 90 mm     |
| Dimensions (W × D)     | 120 mm × 85 mm             |
| Weight                 | 700 g                      |
| Autoclavable           | No                         |
| Pack Size              | 1 Piece (1 Rack × 1 Piece) |
| Color                  | Silver                     |
| Places                 | 8                          |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about Mini Cooler

Small rack with ice-gel tube holder insert for 8 x 1.5/2.0 ml tubes. Perfect for the transportation of a small number of samples.

**NOTE:** Must be frozen to -20 °C for a minimum of 24 hours before use. Not autoclavable.