



Cap-Locks for 1.5 / 2.0 ml Microcentrifuge Tubes



Product Highlights

- Securely seals tubes for rigorous boiling or freezing
- Assorted colours for easier sample identification
- Tab allows tubes to be lifted in/out of boiling water, Dry Baths or freezers
- Writing area allows easy identification or organization into types
- Re-usable
- Compatible with most centrifuge tubes

Product Applications

Seal tubes for rigorous boiling or freezing

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	I1415-1508
Pack Size	100 Pcs. (1 Bag × 100 Pcs.)
For use with (SL)	1.5/2.0 ml Tubes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about Cap-Locks for 1.5 / 2.0 ml Microcentrifuge Tubes

Reusable cap locks to securely seal tubes for rigorous boiling or freezing. Assorted colours to aid sample identification. Assorted Fluorescent: Blue, Green, Red, Orange, Yellow, Black. Compatible with most microcentrifuge tubes.

All Cap-Locks are supplied in mixed packs of eight different colours.