



2.0 ml Screw Cap Tube, Ribbed, Skirted, Natural



Product Highlights

- Optically clear

Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	E1420-2320
Pack Size	500 Tubes (1 Pack × 500 Tubes)
Sterile	No
Autoclavable	Yes
Volume	2.0 ml
Tube Style	Ribbed
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
PCR inhibitor free	Yes
RNase free	Yes
Color	Natural
g-Force max.	20,000 × g
Tube Base	Skirted

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de













More informations about 2.0 ml Screw Cap Tube, Ribbed, Skirted, Natural

Mix and match – tubes and caps!

Polypropylene screw cap tubes that are suitable for many applications from sample storage in the gas phase of liquid nitrogen to freezer, through boiling to autoclave (121 °C).













All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	0.5 ml Screw Cap Tube, Ribbed, Skirted, Natural Volume: 0.5 ml Tube Base: Skirted Tube Style: Ribbed  Create PDF Data Sheet	500 Tubes (1 Pack × 500 Tubes)	E1405-2120
	0.5 ml Screw Cap Tube, Plain, Conical, Natural Volume: 0.5 ml Tube Base: Conical Tube Style: Plain  Create PDF Data Sheet	500 Tubes (1 Pack × 500 Tubes)	E1405-2130
	0.5 ml Screw Cap Tube, Plain, Skirted, Natural Volume: 0.5 ml Tube Base: Skirted Tube Style: Plain  Create PDF Data Sheet	500 Tubes (1 Pack × 500 Tubes)	E1405-2140
	1.5 ml Screw Cap Tube, Ribbed, Conical, Natural Volume: 1.5 ml Tube Base: Conical Tube Style: Ribbed  Create PDF Data Sheet	500 Tubes (1 Pack × 500 Tubes)	E1415-2210
	1.5 ml Screw Cap Tube, Ribbed, Skirted, Natural Volume: 1.5 ml Tube Base: Skirted Tube Style: Ribbed  Create PDF Data Sheet	500 Tubes (1 Pack × 500 Tubes)	E1415-2220

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 1.5 ml Screw Cap Tube, Graduated, Conical, Natural Volume: 1.5 ml Tube Base: Conical Tube Style: Graduated  Create PDF Data Sheet	500 Tubes (1 Pack × 500 Tubes)	E1415-2230
 1.5 ml Screw Cap Tube, Plain, Skirted, Natural Volume: 1.5 ml Tube Base: Skirted Tube Style: Plain  Create PDF Data Sheet	500 Tubes (1 Pack × 500 Tubes)	E1415-2240
 2.0 ml Screw Cap Tube, Ribbed, Skirted, Natural Volume: 2.0 ml Tube Base: Skirted Tube Style: Ribbed  Create PDF Data Sheet	500 Tubes (1 Pack × 500 Tubes)	E1420-2320
 2.0 ml Screw Cap Tube, Graduated, Conical, Natural Volume: 2.0 ml Tube Base: Conical Tube Style: Graduated  Create PDF Data Sheet	500 Tubes (1 Pack × 500 Tubes)	E1420-2330
 2.0 ml Screw Cap Tube, Graduated, Skirted, Natural Volume: 2.0 ml Tube Base: Skirted Tube Style: Graduated  Create PDF Data Sheet	500 Tubes (1 Pack × 500 Tubes)	E1420-2340

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de