



10 µl graduated TipOne® Tip, Natural, Refills (non-sterile)



Product Highlights

- Graduated at 2.5 µl and 10 µl
- Universal tip fit
- High quality virgin polypropylene tips, with extremely smooth surfaces for reduced sample retention
- Certified RNase, DNase, DNA and Pyrogen free
- Refills dramatically reduce plastic waste compared to pre-racked tips
- Refill top and bottom shells made from 100% recycled plastic
- All components, excluding tips and shrink wrap, are recyclable

Product Applications

Pipetting liquids
Dispensing liquids
Mixing liquids
Filling plates and reaction vessels
Phase extraction and removal of supernatants

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	S1111-3700
Sterile	No
Autoclavable	Yes
Filter tip	Non-filter
Volume range	0.1 - 10 µl
Volume	10 µl
Pack Size	960 Tips (10 Refills × 96 Tips)
Color Code Wafer	Red
Tip Type	standard
Material	Polypropylene
Graduations	Yes
Color	Natural
Sustainability attributes	Refills dramatically reduce plastic waste compared to pre-racked tips. Refill top and bottom shells made from 100% recycled plastic. All components, excluding tips and film wrap, are recyclable

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



For use with brand + product (SL)

Abimed: BioPette-E, Discovery Comfort, LABMATE Optima

Biohit: eLine, ePet, mLine, Picus, Proline

BRAND: Transferpette Electronic, Transferpette Pro, Transferpette S

Capp: CappAero Comfort, CappAero Multi, CappAero Universal, CappAero Universal short, ecopipette, Singlechannel

Eppendorf: Reference, Reference 2, Research, Research 3 neo, Research plus, Research pro, Xplorer

Gilson: Pipetman Concept, Pipetman L, Pipetman Neo, Pipetman P, Pipetman Ultra

Hamilton: SoftGrip Pipette, SofTouch

HTL: Discovery Comfort, Discovery Pro, Labmate Pro, Optipette

Nichiryo: 00-NP7, 50F, NLE, NPLO, NPX

Rainin: Classic, EDP1, EDP-3 Plus, Pipet-Lite, Pipet-Plus

Sartorius: eLine, ePet, mLine, Picus, Picus NxT, Proline, Proline Plus, Tacta

Socorex: Acura 815 manual, Acura 825, Acura 855, Acura electro 925, Acura electro 955, Calibra digital

STARLAB: ErgoOne, ErgoOne E, MicroOne, StarPet, StarPet-E

Thermo Fisher Scientific Inc.: Finn Focus, Finn Focus Micro, Finnpipette Digital, Finnpipette F1/F2, Finnpipette Novus, Finnpipette Novus Pink

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about 10 µl graduated TipOne® Tip, Natural, Refills (non-sterile)

High quality, reliable [pipette tips](#) with a universal tip fit – [TipOne](#)® has it all! If you are looking for a sustainable [pipette](#) tip system with accuracy and precision you can trust, [TipOne](#)® is the tip for you.

[TipOne](#)® refills are a quick and easy way to obtain economical racked tips, and they save on plastic use!

Reduce single-use plastic!













Compared to pre-racked tips, using [TipOne](#)® refills dramatically reduces the amount of single-use plastic going through your laboratory. Quick and easy to load, unwrap and click in place into an empty [TipOne](#)® rack. The refills are supplied in two separate towers of 5x96 tips, so you can rack up 5 or 10 at a time. Tip information is printed on every wafer, providing full traceability when racked. Once racked, and if needed, the refills can be autoclaved.

Join our journey to sustainable pipetting

The clear, polypropylene system components are made from 100% recycled plastic. [TipOne](#)® tip racks are built to last and can be reused again and again with these cost-effective refills. All components of [TipOne](#)® refills, excluding tips and shrink wrap, are recyclable. Sign up for our [TipOne](#)® recycling service!














All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 10 µl graduated TipOne® Tip, Natural, Racks (non-sterile) Sterile: No Pack Size: 960 Tips (10 Racks × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Racks × 96 Tips)	S1111-3800
 10 µl graduated TipOne® Tip, Natural, Racks (non-sterile), Case Sterile: No Pack Size: 7,680 Tips (80 Racks × 96 Tips)  Create PDF Data Sheet	7,680 Tips (80 Racks × 96 Tips)	S1111-3800-C
 10 µl graduated TipOne® Tip, Natural, Racks (sterile) Sterile: Yes Pack Size: 960 Tips (10 Racks × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Racks × 96 Tips)	S1111-3810
 10 µl graduated TipOne® Tip, Natural, Racks (sterile), Case Sterile: Yes Pack Size: 7,680 Tips (80 Racks × 96 Tips)  Create PDF Data Sheet	7,680 Tips (80 Racks × 96 Tips)	S1111-3810-C
 10 µl graduated TipOne® Tip, Natural, Refills (non-sterile) Sterile: No Pack Size: 960 Tips (10 Refills × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Refills × 96 Tips)	S1111-3700
 10 µl graduated TipOne® Tip, Natural, Refills (non-sterile), Case Sterile: No Pack Size: 9,600 Tips (100 Refills × 96 Tips)  Create PDF Data Sheet	9,600 Tips (100 Refills × 96 Tips)	S1111-3700-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <p>10 µl graduated TipOne® Tip, Natural, Stack rack (non-sterile) Sterile: No Pack Size: 960 Tips (1 Stack rack × 960 Tips)</p> <p> Create PDF Data Sheet</p>	960 Tips (1 Stack rack × 960 Tips)	S1111-3200
 <p>10 µl graduated TipOne® Tip, Natural, Stack racks (non-sterile), Case Sterile: No Pack Size: 9,600 Tips (10 Stack racks × 960 Tips)</p> <p> Create PDF Data Sheet</p>	9,600 Tips (10 Stack racks × 960 Tips)	S1111-3200-C
 <p>10 µl graduated TipOne® Tip, Natural, Stack rack (sterile) Sterile: Yes Pack Size: 960 Tips (1 Stack rack × 960 Tips)</p> <p> Create PDF Data Sheet</p>	960 Tips (1 Stack rack × 960 Tips)	S1111-3210
 <p>10 µl graduated TipOne® Tip, Natural, Stack racks (sterile), Case Sterile: Yes Pack Size: 9,600 Tips (10 Stack racks × 960 Tips)</p> <p> Create PDF Data Sheet</p>	9,600 Tips (10 Stack racks × 960 Tips)	S1111-3210-C
 <p>10 µl graduated TipOne® Tip, Natural, Bag (non-sterile) Sterile: No Pack Size: 1,000 Tips (1 Bag × 1000 Tips)</p> <p> Create PDF Data Sheet</p>	1,000 Tips (1 Bag × 1000 Tips)	S1111-3000
 <p>10 µl graduated TipOne® Tip, Natural, Bags (non-sterile), Case Sterile: No Pack Size: 50,000 Tips (50 Bags × 1000 Tips)</p>	50,000 Tips (50 Bags × 1000 Tips)	S1111-3000-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
 Neuer Höltigbaum 38
 22143 Hamburg
 Email: info@starlab.de