



## 0.5 ml Screw Cap Tube, Plain, Skirted, Amber



### Product Highlights

- Routinely tested to 20,000 x g
- All tubes and caps are autoclavable
- Certified RNase, DNase and DNA free
- 100% medical grade polypropylene
- Suitable for light sensitive samples

## Product Applications

Suitable for many applications from sample storage in the gas phase of liquid nitrogen to freezer, through boiling to autoclave

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	E1405-2147
Sterile	No
Autoclavable	Yes
Volume	0.5 ml
Pack Size	500 Tubes (1 Bag × 500 Tubes)
Tube Style	Plain
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
PCR inhibitor free	Yes
RNase free	Yes
Color	Amber
g-Force max.	20,000 × g
Tube Base	Skirted

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)













## More informations about 0.5 ml Screw Cap Tube, Plain, Skirted, Amber

### Mix and match – tubes and caps!

Polypropylene screw cap tubes made from 100% medical grade polypropylene. Suitable for many applications from sample storage in the gas phase of liquid nitrogen to freezer, through boiling to autoclave (121 °C). Available with or without graduations, the graduated versions are marked at every 100 µl to provide easy sample estimation. Offered with either skirted or conical bases. Routinely tested to 20,000 x g. Caps are available separately; all tubes and caps are autoclavable. Certified RNase, DNase and DNA free.







## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">0.5 ml Screw Cap Tube, Ribbed, Skirted, Amber</a> <b>Volume:</b> 0.5 ml <b>Tube Style:</b> Ribbed <b>Tube Base:</b> Skirted   <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Bag × 500 Tubes)	E1405-2127
 <a href="#">1.5 ml Screw Cap Tube, Ribbed, Skirted, Amber</a> <b>Volume:</b> 1.5 ml <b>Tube Style:</b> Ribbed <b>Tube Base:</b> Skirted   <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Bag × 500 Tubes)	E1415-2227
 <a href="#">1.5 ml Screw Cap Tube, Graduated, Conical, Amber</a> <b>Volume:</b> 1.5 ml <b>Tube Style:</b> Graduated <b>Tube Base:</b> Conical   <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Bag × 500 Tubes)	E1415-2237
 <a href="#">2.0 ml Screw Cap Tube, Graduated, Conical, Amber</a> <b>Volume:</b> 2.0 ml <b>Tube Style:</b> Graduated <b>Tube Base:</b> Conical   <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Bag × 500 Tubes)	E1420-2337
 <a href="#">2.0 ml Screw Cap Tube, Graduated, Skirted, Amber</a> <b>Volume:</b> 2.0 ml <b>Tube Style:</b> Graduated <b>Tube Base:</b> Skirted   <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Bag × 500 Tubes)	E1420-2347

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">1.5 ml Screw Cap Tube, Plain, Skirted, Amber</a> <b>Volume:</b> 1.5 ml <b>Tube Style:</b> Plain <b>Tube Base:</b> Skirted   <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Bag × 500 Tubes)	E1415-2247
 <a href="#">0.5 ml Screw Cap Tube, Plain, Skirted, Amber</a> <b>Volume:</b> 0.5 ml <b>Tube Style:</b> Plain <b>Tube Base:</b> Skirted   <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Bag × 500 Tubes)	E1405-2147

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)