



## 0.2 ml 8-Strip PCR Tubes, Individually Attached Domed Caps, Natural



### Product Highlights

- Manufactured from medical grade polypropylene
- Attached caps eliminates the possibility of cross-contamination through accidental tube–cap mix up
- PCR inhibitor free
- Certified RNase, Dnase, DNA and Pyrogen-free
- Ultra thin walls ensure efficient thermal transfer and maximum yield
- Supplied in resealable bags

## Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

|                    |                                                                 |
|--------------------|-----------------------------------------------------------------|
| Art. No.           | I1402-2900                                                      |
| Pack Size          | 120 Tube Strips with Caps (1 Packs × 120 Tube Strips with Caps) |
| Sterile            | No                                                              |
| Autoclavable       | Yes                                                             |
| Volume             | 0.2 ml                                                          |
| No. of tubes       | 8                                                               |
| Cap style          | Individually attached domed caps                                |
| Material           | Polypropylene                                                   |
| DNA free           | Yes                                                             |
| DNase free         | Yes                                                             |
| Free of endotoxins | Yes                                                             |
| Human DNA free     | Yes                                                             |
| PCR inhibitor free | Yes                                                             |
| Pyrogen free       | Yes                                                             |
| RNase free         | Yes                                                             |
| Color              | Natural                                                         |
| Suitable for PCR   | yes                                                             |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about 0.2 ml 8-Strip PCR Tubes, Individually Attached Domed Caps, Natural

Starlab PCR Strips are manufactured from medical grade polypropylene. The thin-walled tubes ensure efficient and homogenous heat transfer. The attached 'crocodile' style strip caps eliminate the possibility of cross-contamination from cap mix ups. Available in 8- or 12-Strip Tubes and Caps.

Low-Profile Strips (0.1 ml): 15.24 mm high, for use in cyclers with low-profile capability only.

Standard Strips (0.2 ml): 20.83 mm high, for use in the majority of cyclers.



## All Variations

| PRODUCT NAME                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       | PACKAGING SIZE                                                        | ART. NO.     |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------|--------------|
| <br><a href="#">0.2 ml 8-Strip PCR Tubes, Individually Attached Flat Caps (Xtra-Clear)</a><br><b>Cap style:</b> Individually attached flat caps (X-Clear)<br><b>Volume:</b> 0.2 ml<br><b>Pack Size:</b> 120 Tube Strips with Caps (1 Pack × 120 Tube Strips with Caps)<br> <a href="#">Create PDF Data Sheet</a>                 | 120 Tube Strips with Caps<br>(1 Pack × 120 Tube Strips with Caps)     | A1402-3700   |
| <br><a href="#">0.2 ml 8-Strip PCR Tubes, Individually Attached Flat Caps (Xtra-Clear)(Case Size)</a><br><b>Cap style:</b> Individually attached flat caps (X-Clear)<br><b>Volume:</b> 0.2 ml<br><b>Pack Size:</b> 1,200 Tube Strips with Caps (10 Packs × 120 Tube Strips with Caps)<br> <a href="#">Create PDF Data Sheet</a> | 1,200 Tube Strips with Caps<br>(10 Packs × 120 Tube Strips with Caps) | A1402-3700-C |
| <br><a href="#">0.2 ml 8-Strip PCR Tubes, Individually Attached Domed Caps, Natural</a><br><b>Cap style:</b> Individually attached domed caps<br><b>Volume:</b> 0.2 ml<br><b>Pack Size:</b> 120 Tube Strips with Caps (1 Packs × 120 Tube Strips with Caps)<br> <a href="#">Create PDF Data Sheet</a>                        | 120 Tube Strips with Caps<br>(1 Packs × 120 Tube Strips with Caps)    | I1402-2900   |
| <br><a href="#">8-Strip PCR Tubes with Individual Attached Caps, 0.2 ml (mixed)</a><br><b>Cap style:</b> Individually attached domed caps<br><b>Volume:</b> 0.2 ml<br><b>Pack Size:</b> 120 Tube Strips with Caps (5 Packs × 24 Tube Strips with Caps)<br> <a href="#">Create PDF Data Sheet</a>                             | 120 Tube Strips with Caps<br>(5 Packs × 24 Tube Strips with Caps)     | I1402-2908   |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



| PRODUCT NAME                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      | PACKAGING SIZE                                                        | ART. NO.     |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------|--------------|
| <br><a href="#">0.1 ml 8-Strip PCR Tubes, Low-Profile, Individually Attached Flat Caps (Xtra-Clear)</a><br><b>Cap style:</b> Individually attached flat caps – (Xtra-Clear), low-profile<br><b>Volume:</b> 0.1 ml<br><b>Pack Size:</b> 120 Tube Strips with Caps (1 Pack × 120 Tube Strips with Caps)<br> <a href="#">Create PDF Data Sheet</a> | 120 Tube Strips with Caps<br>(1 Pack × 120 Tube Strips with Caps)     | A1402-3800   |
| <br><a href="#">0.2 ml 8-Strip PCR Tubes, Individually Attached Domed Caps, Natural (Case Size)</a><br><b>Cap style:</b> Individually attached domed caps<br><b>Volume:</b> 0.2 ml<br><b>Pack Size:</b> 1,200 Tube Strips with Caps (10 Packs × 120 Tube Strips with Caps)<br> <a href="#">Create PDF Data Sheet</a>                            | 1,200 Tube Strips with Caps<br>(10 Packs × 120 Tube Strips with Caps) | I1402-2900-C |
| <br><a href="#">0.2ml 8-Strip PCR Tubes, Individually Attached Domed Shielded Caps</a><br><b>Cap style:</b> Individually attached domed shielded caps<br><b>Volume:</b> 0.2 ml<br><b>Pack Size:</b> 120 Tube Strips with Caps (1 Pack × 120 Tube Strips with Caps)<br> <a href="#">Create PDF Data Sheet</a>                                | 120 Tube Strips with Caps<br>(1 Pack × 120 Tube Strips with Caps)     | S1602-2900   |
| <br><a href="#">0.2 ml 8-Strip PCR Tubes, Natural with Attached Flat (Xtra-Clear) Cap Strip</a><br><b>Cap style:</b> Attached flat cap strip (Xtra-Clear)<br><b>Volume:</b> 0.2 ml<br><b>Pack Size:</b> 125 Tube Strips with Caps (1 Pack × 125 Tube Strips with Caps)<br> <a href="#">Create PDF Data Sheet</a>                            | 125 Tube Strips with Caps<br>(1 Pack × 125 Tube Strips with Caps)     | I1402-2700   |
| <br><a href="#">0.2 ml 8-Strip PCR Tubes, Natural with Attached Domed Cap Strip</a><br><b>Cap style:</b> Attached domed cap strip<br><b>Volume:</b> 0.2 ml<br><b>Pack Size:</b> 125 Tube Strips with Caps (1 Pack × 125 Tube Strips with Caps)<br> <a href="#">Create PDF Data Sheet</a>                                                    | 125 Tube Strips with Caps<br>(1 Pack × 125 Tube Strips with Caps)     | I1402-2800   |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



| PRODUCT NAME                                                                      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                               | PACKAGING SIZE                                                        | ART. NO.     |
|-----------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------|--------------|
|  | <a href="#">0.2 ml 8-Strip PCR Tubes, Natural with Attached Flat (Xtra-Clear) Cap Strip (Case Size)</a><br><b>Cap style:</b> Attached flat cap strip (Xtra-Clear)<br><b>Volume:</b> 0.2 ml<br><b>Pack Size:</b> 1,250 Tube Strips with Caps<br>(10 Packs × 125 Tube Strips with Caps)<br> <a href="#">Create PDF Data Sheet</a>                                                              | 1,250 Tube Strips with Caps<br>(10 Packs × 125 Tube Strips with Caps) | I1402-2700-C |
|                                                                                   | <br><a href="#">0.2 ml 8-Strip PCR Tubes, Natural with Attached Domed Cap Strip (Case Size)</a><br><b>Cap style:</b> Attached domed cap strip<br><b>Volume:</b> 0.2 ml<br><b>Pack Size:</b> 1,250 Tube Strips with Caps<br>(10 Packs × 125 Tube Strips with Caps)<br> <a href="#">Create PDF Data Sheet</a> |                                                                       |              |