



0.2 ml 8-Strip PCR Tubes, Natural with Attached Domed Cap Strip



Product Highlights

- Manufactured from medical grade polypropylene
- Individually attached caps eliminates cross-contamination from cap mix-up
- PCR inhibitor free
- Certified RNase, Dnase, DNA and Pyrogen-free
- Ultra thin walls ensure efficient thermal transfer and maximum yield
- Supplied in resealable bags

Product Applications

Standard PCR (Endpoint PCR)

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	I1402-2800
Sterile	No
Autoclavable	Yes
Volume	0.2 ml
Pack Size	125 Tube Strips with Caps (1 Pack × 125 Tube Strips with Caps)
No. of tubes	8
Cap style	Attached domed cap strip
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
Human DNA free	Yes
PCR inhibitor free	Yes
Pyrogen free	Yes
RNase free	Yes
Color	Natural
Suitable for PCR	yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about 0.2 ml 8-Strip PCR Tubes, Natural with Attached Domed Cap Strip









Help eliminate cross-contamination with [PCR Strips](#) with Attached Cap Strip

Discover convenience and efficiency with Starlab's [PCR strips](#) with attached cap strip, designed for mid-throughput settings. These [PCR strips](#) ensure all tubes are equally protected, supporting consistent results across samples, as well as simplifying handling to increase workflow efficiency. The attached 'crocodile' style strip caps eliminate the possibility of cross-contamination from cap mix ups.

The 0.2ml tubes with attached domed cap strip are suitable for use in the majority of cyclers. Their ultra-thin, uniform walls ensure optimal heat transfer. Certified free of RNase, DNase, DNA and PCR inhibitors, and tested Pyrogen free.













All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 0.2 ml 8-Strip PCR Tubes, Individually Attached Flat Caps (Xtra-Clear) Cap style: Individually attached flat caps (X-Clear) Volume: 0.2 ml Pack Size: 120 Tube Strips with Caps (1 Pack × 120 Tube Strips with Caps)  Create PDF Data Sheet	120 Tube Strips with Caps (1 Pack × 120 Tube Strips with Caps)	A1402-3700
 0.2 ml 8-Strip PCR Tubes, Individually Attached Flat Caps (Xtra-Clear)(Case Size) Cap style: Individually attached flat caps (X-Clear) Volume: 0.2 ml Pack Size: 1,200 Tube Strips with Caps (10 Packs × 120 Tube Strips with Caps)  Create PDF Data Sheet	1,200 Tube Strips with Caps (10 Packs × 120 Tube Strips with Caps)	A1402-3700-C
 0.2 ml 8-Strip PCR Tubes, Individually Attached Domed Caps, Natural Cap style: Individually attached domed caps Volume: 0.2 ml Pack Size: 120 Tube Strips with Caps (1 Packs × 120 Tube Strips with Caps)  Create PDF Data Sheet	120 Tube Strips with Caps (1 Packs × 120 Tube Strips with Caps)	I1402-2900
 8-Strip PCR Tubes with Individual Attached Caps, 0.2 ml (mixed) Cap style: Individually attached domed caps Volume: 0.2 ml Pack Size: 120 Tube Strips with Caps (5 Packs × 24 Tube Strips with Caps)  Create PDF Data Sheet	120 Tube Strips with Caps (5 Packs × 24 Tube Strips with Caps)	I1402-2908

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de







PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <p>0.1 ml 8-Strip PCR Tubes, Low-Profile, Individually Attached Flat Caps (Xtra-Clear) Cap style: Individually attached flat caps – (Xtra-Clear), low-profile Volume: 0.1 ml Pack Size: 120 Tube Strips with Caps (1 Pack × 120 Tube Strips with Caps)</p> <p> Create PDF Data Sheet</p>	<p>120 Tube Strips with Caps (1 Pack × 120 Tube Strips with Caps)</p>	<p>A1402-3800</p>
 <p>0.2 ml 8-Strip PCR Tubes, Individually Attached Domed Caps, Natural (Case Size) Cap style: Individually attached domed caps Volume: 0.2 ml Pack Size: 1,200 Tube Strips with Caps (10 Packs × 120 Tube Strips with Caps)</p> <p> Create PDF Data Sheet</p>	<p>1,200 Tube Strips with Caps (10 Packs × 120 Tube Strips with Caps)</p>	<p>I1402-2900-C</p>
 <p>0.2ml 8-Strip PCR Tubes, Individually Attached Domed Shielded Caps Cap style: Individually attached domed shielded caps Volume: 0.2 ml Pack Size: 120 Tube Strips with Caps (1 Pack × 120 Tube Strips with Caps)</p> <p> Create PDF Data Sheet</p>	<p>120 Tube Strips with Caps (1 Pack × 120 Tube Strips with Caps)</p>	<p>S1602-2900</p>
 <p>0.2 ml 8-Strip PCR Tubes, Natural with Attached Flat (Xtra-Clear) Cap Strip Cap style: Attached flat cap strip (Xtra-Clear) Volume: 0.2 ml Pack Size: 125 Tube Strips with Caps (1 Pack × 125 Tube Strips with Caps)</p> <p> Create PDF Data Sheet</p>	<p>125 Tube Strips with Caps (1 Pack × 125 Tube Strips with Caps)</p>	<p>I1402-2700</p>
 <p>0.2 ml 8-Strip PCR Tubes, Natural with Attached Domed Cap Strip Cap style: Attached domed cap strip Volume: 0.2 ml Pack Size: 125 Tube Strips with Caps (1 Pack × 125 Tube Strips with Caps)</p> <p> Create PDF Data Sheet</p>	<p>125 Tube Strips with Caps (1 Pack × 125 Tube Strips with Caps)</p>	<p>I1402-2800</p>

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
 Neuer Höltigbaum 38
 22143 Hamburg
 Email: info@starlab.de



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <p>0.2 ml 8-Strip PCR Tubes, Natural with Attached Flat (Xtra-Clear) Cap Strip (Case Size) Cap style: Attached flat cap strip (Xtra-Clear) Volume: 0.2 ml Pack Size: 1,250 Tube Strips with Caps (10 Packs × 125 Tube Strips with Caps)</p> <p> Create PDF Data Sheet</p>	1,250 Tube Strips with Caps (10 Packs × 125 Tube Strips with Caps)	I1402-2700-C
 <p>0.2 ml 8-Strip PCR Tubes, Natural with Attached Domed Cap Strip (Case Size) Cap style: Attached domed cap strip Volume: 0.2 ml Pack Size: 1,250 Tube Strips with Caps (10 Packs × 125 Tube Strips with Caps)</p> <p> Create PDF Data Sheet</p>	1,250 Tube Strips with Caps (10 Packs × 125 Tube Strips with Caps)	I1402-2800-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de