



0.2 ml 8-Strip PCR Tubes, White with Attached Flat (Xtra-Clear) Cap Strip



Product Highlights

- Manufactured from medical grade polypropylene
- Individually attached caps eliminates cross-contamination from cap mix-up
- PCR inhibitor free
- Suitable for Realtime PCR (quantitative PCR = qPCR)
- Certified RNase, Dnase, DNA and Pyrogen-free
- Ultra thin walls ensure efficient thermal transfer and maximum yield
- Supplied in resealable bags

Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	I1409-2700
Pack Size	125 Tube Strips with Caps (1 Pack × 125 Tube Strips with Caps)
Sterile	No
Autoclavable	Yes
Volume	0.2 ml
No. of tubes	8
Cap style	Attached flat cap strip (Xtra-Clear)
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
Human DNA free	Yes
PCR inhibitor free	Yes
Pyrogen free	Yes
RNase free	Yes
Color	White
Suitable for PCR	yes
Suitable for Real-Time PCR (qPCR)	yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about 0.2 ml 8-Strip PCR Tubes, White with Attached Flat (Xtra-Clear) Cap Strip

Starlab PCR Strips are manufactured from medical grade polypropylene. The thin-walled tubes ensure efficient and homogenous heat transfer. The attached 'crocodile' style strip caps eliminate the possibility of cross-contamination from cap mix ups. Available in 8- or 12-Strip Tubes and Caps.

Low-Profile Strips (0.1 ml): 15.24 mm high, for use in cyclers with low-profile capability only.

Standard Strips (0.2 ml): 20.83 mm high, for use in the majority of cyclers.