



## 96-Well PCR Plate, Semi-Skirted



### Product Highlights

- Approximate maximum capacity for all 96-well plates is 350 µl for standard height plates, or 200 µl for low-profile plates
- Available as natural, white or black plates. White plates are ideal for qPCR, giving optimal signal-to-noise ratio for fluorescent based assays. Black plates help to minimize light diffusion.
- Specialist plates (for FAST® Systems, qPCR)

## Product Applications

PCR

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	I1402-9700
Material	Polypropylene
Profile	Standard height
Elevated wells	no
Raised rim	no
Matrix	Printed
Pack Size	10 Plates (1 Box × 10 Plates)
Free of heavy metals	Yes
Free of endotoxins	Yes
Pyrogen free	Yes
Cut corner	A12
DNA free	Yes
PCR inhibitor free	Yes
RNase free	Yes
DNase free	Yes
Suitable for Real-Time PCR (qPCR)	yes
Suitable for PCR	yes
Autoclavable	Yes
Sterile	No
Number of wells	96 wells
Frame Design	Semi-skirted
Color	Natural

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



Volume

350 µl

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about 96-Well PCR Plate, Semi-Skirted


Starlab Semi-Skirted PCR Plates feature a 7 mm high skirt around the complete plate edge.

The skirt improves the stability of the plate and makes them more suited e. g. for automated applications. Furthermore a labelling on the plate edge e. g. by barcode-labels is possible. Plates with elevated wells reduce the contamination risk by cross-contamination.





- › Autoclavable
- › Supplied in resealable bags
- › Convenient small pack sizes
- › Certified RNase, DNase, DNA and Pyrogen free
- › PCR Inhibitor free

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## Accessories

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	<a href="#">Polypropylene PCR Sealing Film Strips, Clear</a>	200 Pcs. (1 Box × 200 Pcs.)	E2796-2850
	<a href="#">Polyester PCR Sealing Film, Clear</a>	100 Pcs. (1 Box × 100 Pcs.)	E2796-0100
	<a href="#">Xtra-Clear Advanced Polyolefin StarSeal (qPCR)</a>	100 Pcs. (1 Box × 100 Pcs.)	E2796-9795
	<a href="#">Clear Polypropylene Seal (PCR)</a>	100 Pcs. (1 Box × 100 Pcs.)	E2796-0793
	<a href="#">Clear Polyolefin StarSeal (PCR)</a>	100 Pcs. (1 Box × 100 Pcs.)	E2796-9793
	<a href="#">Aluminium Sealing Film, 60 µm (PCR)</a>	100 Pcs. (1 Box × 100 Pcs.)	E2796-0792
	<a href="#">Aluminium StarSeal (PCR)</a>	100 Pcs. (1 Box × 100 Pcs.)	E2796-9792
	<a href="#">Silicone Sealing Mat for 96-Well PCR Plate</a>	5 Pcs. (1 Box × 5 Pcs.)	E1403-0000

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)