



384-Well Skirted PCR Plates for Roche® Lightcycler®, White



Product Highlights

- Approximate maximum capacity for all 96-well plates is 350 µl for standard height plates, or 200 µl for low-profile plates
- Holes on sides for robotics handling
- Certified RNase, DNase, DNA and Pyrogen-free
- Opaque white plates available for qPCR
- PCR Inhibitor free
- Supplied in resealable bags

Product Applications

PCR

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

| | |
|-----------------------------------|--------------------------------|
| Art. No. | E1042-9909 |
| Sterile | No |
| Autoclavable | Yes |
| Volume | 35 µl |
| Pack Size | 10 Plates (1 Pack × 10 Plates) |
| Number of wells | 384 wells |
| Material | Polypropylene |
| DNA free | Yes |
| DNase free | Yes |
| Free of endotoxins | Yes |
| PCR inhibitor free | Yes |
| Pyrogen free | Yes |
| RNase free | Yes |
| Color | White |
| Cut corner | A24, P24 |
| Frame Design | Skirted |
| Matrix | Printed |
| Profile | For Roche® Lightcycler |
| Suitable for PCR | yes |
| Suitable for Real-Time PCR (qPCR) | yes |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about 384-Well Skirted PCR Plates for Roche® Lightcycler®, White








Starlab Skirted 96- and 384-Well PCR Plates feature a 15 mm high full skirt around the complete plate edge. The skirt provides an extremely high plate stability, which is beneficial specifically in automated "High-Throughput" applications. These type of plates offer the largest area for marking on edge of the plates e. g. by barcodes. Starlab Skirted PCR Plates are available in standard or low profile format.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



Accessories

| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|--------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------|------------|
|  Polyester PCR Sealing Film, Clear | 100 Pcs. (1 Box × 100 Pcs.) | E2796-0100 |
|  Polypropylene PCR Sealing Film Strips, Clear | 200 Pcs. (1 Box × 200 Pcs.) | E2796-2850 |
|  Xtra-Clear Advanced Polyolefin StarSeal (qPCR) | 100 Pcs. (1 Box × 100 Pcs.) | E2796-9795 |
|  Clear Polypropylene Seal (PCR) | 100 Pcs. (1 Box × 100 Pcs.) | E2796-0793 |
|  Clear Polyolefin StarSeal (PCR) | 100 Pcs. (1 Box × 100 Pcs.) | E2796-9793 |
|  Aluminium Sealing Film, 60 µm (PCR) | 100 Pcs. (1 Box × 100 Pcs.) | E2796-0792 |
|  Aluminium StarSeal (PCR) | 100 Pcs. (1 Box × 100 Pcs.) | E2796-9792 |
|  Silicone Sealing Mat for 96-Well PCR Plate | 5 Pcs. (1 Box × 5 Pcs.) | E1403-0000 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de