



384-Well PCR Plate for Roche Lightcycler, Skirted, with Barcode



Product Highlights

- Code 128 Barcoding plus human-readable number
- All barcodes seven digits long
- Low, minimum ordering quantity of 100 plates
- Certified RNase, DNase, DNA and Pyrogen-free
- PCR Inhibitor free
- 96-well and 384-well plate options

Product Applications

PCR

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	E1042-9900-BC
Sterile	No
Autoclavable	Yes
Volume	35 µl
Pack Size	100 Plates (10 Packs × 10 Plates)
Number of wells	384 wells
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
PCR inhibitor free	Yes
Pyrogen free	Yes
RNase free	Yes
Color	Natural
Cut corner	A24, P24
Elevated wells	no
Frame Design	Skirted
Profile	For Roche® Lightcycler
Raised rim	no
Suitable for PCR	yes
Suitable for Real-Time PCR (qPCR)	yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about 384-Well PCR Plate for Roche Lightcycler, Skirted, with Barcode

Customised barcoding is also available












If you require a specific sequence start to your barcoding or require barcoding on plates other than those shown above, we can provide you with a quote for customised barcoding. Note that the minimum order quantity for custom barcoding is 2000 plates. Contact Starlab or your Starlab Account Manager for more details.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 384-Well PCR Plate, ABI/Universal-Style, Skirted, with Barcode Profile: ABI®/Universal-Style Frame Design: Skirted	100 Plates (10 Packs × 10 Plates)	E1042-3840-BC
 384-Well PCR Plate for Roche Lightcycler, Skirted, with Barcode Profile: For Roche® Lightcycler Frame Design: Skirted  Create PDF Data Sheet	100 Plates (10 Packs × 10 Plates)	E1042-9900-BC
 96-Well PCR Plate, Semi-Skirted, Raised Rim, with Barcode Profile: Raised rim Frame Design: Semi-skirted  Create PDF Data Sheet	100 Plates (10 Packs × 10 Plates)	E1403-8200-BC
 96-Well PCR Plate, White, Semi-Skirted, for Roche® Lightcycler, with Barcode Profile: Low profile, for Roche® Lightcycler® Frame Design: Semi-skirted  Create PDF Data Sheet	100 Plates (10 Packs × 10 Plates)	I1402-9909-BC
 96-Well PCR Plate, Skirted, Low-Profile, with Barcode Profile: Low Frame Design: Skirted  Create PDF Data Sheet	100 Plates (10 Packs × 10 Plates)	E1403-5200-BC
 96-Well PCR Plate, Semi-Skirted, Straight Edges, with Barcode Profile: Standard height, straight edges Frame Design: Semi-skirted  Create PDF Data Sheet	100 Plates (10 Packs × 10 Plates)	I1402-9800-BC

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



Accessories

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 Silicone Sealing Mat for 96-Well PCR Plate	5 Pcs. (1 Box × 5 Pcs.)	E1403-0000
 Aluminium StarSeal (PCR)	100 Pcs. (1 Box × 100 Pcs.)	E2796-9792
 Aluminium Sealing Film, 60 µm (PCR)	100 Pcs. (1 Box × 100 Pcs.)	E2796-0792
 Clear Polyolefin StarSeal (PCR)	100 Pcs. (1 Box × 100 Pcs.)	E2796-9793
 Clear Polypropylene Seal (PCR)	100 Pcs. (1 Box × 100 Pcs.)	E2796-0793
 Xtra-Clear Advanced Polyolefin StarSeal (qPCR)	100 Pcs. (1 Box × 100 Pcs.)	E2796-9795
 Polyester PCR Sealing Film, Clear	100 Pcs. (1 Box × 100 Pcs.)	E2796-0100
 Polypropylene PCR Sealing Film Strips, Clear	200 Pcs. (1 Box × 200 Pcs.)	E2796-2850

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de