



0.1 ml 8-Strip Non-Flex PCR Tubes, Low-Profile, White, Individually Attached Flat Caps (Xtra-Clear)



Product Highlights

- Manufactured from medical grade polypropylene
- 'Non-flex' PCR Strips for easier handling
- Xtra-Clear caps for qPCR
- Certified RNase, Dnase, DNA & Pyrogen-free
- PCR Inhibitor free
- Ultra thin walls ensure efficient thermal transfer and maximum yield
- Supplied in resealable bags
- Opaque white tubes to maximise reflectance for enhanced signal

Product Applications

Realtime PCR (quantitative PCR = qPCR)

–

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	I1409-3800
Sterile	No
Autoclavable	Yes
Volume	0.1 ml
Pack Size	120 Tube Strips with Caps (1 Bag × 120 Tube Strips with Caps)
Cap style	Individually attached flat caps – (Xtra-Clear), low-profile
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
Human DNA free	Yes
PCR inhibitor free	Yes
Pyrogen free	Yes
RNase free	Yes
Color	White
Profile	Low
Suitable for PCR	yes
Suitable for Real-Time PCR (qPCR)	yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about 0.1 ml 8-Strip Non-Flex PCR Tubes, Low-Profile, White, Individually Attached Flat Caps (Xtra-Clear)











Starlab non-flex tubes have a robust, triple link between each tube making them less flexible and easier to handle. The individually attached caps eliminates possible cross contamination. The caps are set an angle which allows the tubes to be nested together in racks without the lids interfering with one another.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 0.2 ml 8-Strip Non-Flex PCR Tubes, Natural, Individually Attached Flat Caps (Xtra-Clear) Volume: 0.2 ml Color: Natural  Create PDF Data Sheet	120 Tube Strips with Caps (1 Bag × 120 Tube Strips with Caps)	I1402-3700
 0.1 ml 8-Strip Non-Flex PCR Tubes, Low-Profile, Natural, Individually Attached Flat Caps (Xtra-Clear) Volume: 0.1 ml Color: Natural  Create PDF Data Sheet	120 Tube Strips with Caps (1 Bag × 120 Tube Strips with Caps)	I1402-3800
 0.2 ml 8-Strip Non-Flex PCR Tubes, White, Individually Attached Flat Caps (Xtra-Clear) Volume: 0.2 ml Color: White  Create PDF Data Sheet	120 Tube Strips with Caps (1 Bag × 120 Tube Strips with Caps)	I1409-3700
 0.1 ml 8-Strip Non-Flex PCR Tubes, Low-Profile, White, Individually Attached Flat Caps (Xtra-Clear) Volume: 0.1 ml Color: White  Create PDF Data Sheet	120 Tube Strips with Caps (1 Bag × 120 Tube Strips with Caps)	I1409-3800
 0.2 ml 8-Strip Non-Flex PCR Tubes, Mixed, Individually Attached Flat Caps (Xtra-Clear) Volume: 0.2 ml Color: Mixed  Create PDF Data Sheet	120 Tube Strips with Caps (1 Bag × 120 Tube Strips with Caps)	I1402-3708

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de