



## 0.2 ml 8-Strip Non-Flex PCR Tubes, Mixed, Individually Attached Flat Caps (Xtra-Clear)



### Product Highlights

- Manufactured from medical grade polypropylene
- 'Non-flex' PCR Strips for easier handling
- Xtra-Clear caps for qPCR
- Certified RNase, Dnase, DNA & Pyrogen-free
- PCR Inhibitor free
- Ultra thin walls ensure efficient thermal transfer and maximum yield
- Supplied in resealable bags
- Mixed colours for sample classification

## Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

|                                   |   |
|-----------------------------------|---|
| Art. No.                          | I1402-3708  |
| Pack Size                         | 120 Tube Strips with Caps (1 Bag × 120 Tube Strips with Caps) |
| Sterile                           | No  |
| Autoclavable                      | Yes   |
| Volume                            | 0.2 ml  |
| Cap style                         | Individually attached flat caps (Xtra-Clear)                  |
| Material                          | Polypropylene   |
| DNA free                          | Yes   |
| DNase free                        | Yes   |
| Free of endotoxins                | Yes   |
| Human DNA free                    | Yes   |
| PCR inhibitor free                | Yes   |
| Pyrogen free                      | Yes   |
| RNase free                        | Yes   |
| Color                             | Mixed   |
| Profile                           | Standard height   |
| Suitable for PCR                  | yes   |
| Suitable for Real-Time PCR (qPCR) | yes   |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)













## More informations about 0.2 ml 8-Strip Non-Flex PCR Tubes, Mixed, Individually Attached Flat Caps (Xtra-Clear)

Starlab non-flex tubes have a robust, triple link between each tube making them less flexible and easier to handle. The individually attached caps eliminates possible cross contamination. The caps are set an angle which allows the tubes to be nested together in racks without the lids interfering with one another.



## All Variations

| PRODUCT NAME  | PACKAGING SIZE  | ART. NO.   |
|---|---|------------|
| <br><a href="#">0.2 ml 8-Strip Non-Flex PCR Tubes, Natural, Individually Attached Flat Caps (Xtra-Clear)</a><br><b>Volume:</b> 0.2 ml<br><b>Color:</b> Natural<br> <a href="#">Create PDF Data Sheet</a>              | 120 Tube Strips with Caps (1 Bag × 120 Tube Strips with Caps) | I1402-3700 |
| <br><a href="#">0.1 ml 8-Strip Non-Flex PCR Tubes, Low-Profile, Natural, Individually Attached Flat Caps (Xtra-Clear)</a><br><b>Volume:</b> 0.1 ml<br><b>Color:</b> Natural<br> <a href="#">Create PDF Data Sheet</a> | 120 Tube Strips with Caps (1 Bag × 120 Tube Strips with Caps) | I1402-3800 |
| <br><a href="#">0.2 ml 8-Strip Non-Flex PCR Tubes, White, Individually Attached Flat Caps (Xtra-Clear)</a><br><b>Volume:</b> 0.2 ml<br><b>Color:</b> White<br> <a href="#">Create PDF Data Sheet</a>               | 120 Tube Strips with Caps (1 Bag × 120 Tube Strips with Caps) | I1409-3700 |
| <br><a href="#">0.1 ml 8-Strip Non-Flex PCR Tubes, Low-Profile, White, Individually Attached Flat Caps (Xtra-Clear)</a><br><b>Volume:</b> 0.1 ml<br><b>Color:</b> White<br> <a href="#">Create PDF Data Sheet</a> | 120 Tube Strips with Caps (1 Bag × 120 Tube Strips with Caps) | I1409-3800 |
| <br><a href="#">0.2 ml 8-Strip Non-Flex PCR Tubes, Mixed, Individually Attached Flat Caps (Xtra-Clear)</a><br><b>Volume:</b> 0.2 ml<br><b>Color:</b> Mixed<br> <a href="#">Create PDF Data Sheet</a>              | 120 Tube Strips with Caps (1 Bag × 120 Tube Strips with Caps) | I1402-3708 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)