



Sealing Mat for 96-Well Plates with Square Wells, Piercable



Product Highlights

- The chemically resistant mats can be used with DMSO and other solvents.
- Autoclavable
- Piercable

Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	E2896-1801
Pack Size	10 Pcs. (1 Box × 10 Pcs.)
Number of wells	96 wells
Color	Natural
Max. temp.	120 °C
Temperature range	-80 to 120 °C
Temperature min. (SL)	-80 °C
Resistant to DMSO	no
Suitable for Real-Time PCR (qPCR)	no

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about Sealing Mat for 96-Well Plates with Square Wells, Piercable

Particularly suitable for sealing plates with defined heated-lid contact pressure, silicone mats are very flexible, autoclavable and reusable. The mats can also be cut to meet individual needs. Sealing mats are not suitable for Real-Time applications. For best results, we recommend using only the dedicated mat for your plate.