




StarTub PVC Reagent Reservoir for Multichannel Pipettes



Product Highlights

- Clear PVC Reservoir
- V-shaped bottom
- For single use
- Capacity: 55 ml

Product Applications

Reagent takeup with multichannel pipette

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

| | |
|--------------|-----------------------------|
| Art. No. | E2310-1000 |
| Pack Size | 100 Pcs. (1 Box × 100 Pcs.) |
| Autoclavable | No |
| Volume | 55 ml |
| Bottoms | Conical |
| Reusable | No |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about StarTub PVC Reagent Reservoir for Multichannel Pipettes

Disposable, clear PVC reservoirs with V-shaped bottom for single use.