




## 12-Channel Reservoir



### Product Highlights

- Each trough is indexed with a moulded number
- Autoclavable
- Certified RNase, DNase, DNA and PCR Inhibitor free

## Product Applications

Reagent takeup with 12-channel pipette

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

|                        |                          |
|------------------------|--------------------------|
| Art. No.               | E2310-1200               |
| Dimensions (W × D × H) | 128 mm × 87 mm × 44.5 mm |
| Dimensions (W × D)     | 128 mm × 87 mm           |
| Autoclavable           | Yes                      |
| Volume                 | 264 ml                   |
| Pack Size              | 5 Pcs. (1 Box × 5 Pcs.)  |
| Number of troughs      | 12                       |
| Material               | Polypropylene            |
| DNA free               | Yes                      |
| DNase free             | Yes                      |
| PCR inhibitor free     | Yes                      |
| RNase free             | Yes                      |
| Color                  | Clear                    |
| Bottoms                | Pyramid                  |
| Reusable               | Yes                      |
| For use with (SL)      | 12-Channel Pipette       |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about 12-Channel Reservoir

Reusable, clear polypropylene reservoir with 12 x 22 ml troughs, each with 8 x pyramid bottoms.

Dimensions (WxDxH): 128 x 87 x 44.5 mm

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)