




3.0 ml Fine Tip Pasteur Pipette (Sterile), Ind. Wrapped



Product Highlights

- Low affinity surface
- Attached bulb avoids cross-contamination
- One-step biohazard disposal
- Shatter-proof low-density polyethylene, which is non-toxic and inert
- Can be frozen in liquid nitrogen
- Sterile pipettes supplied in peel-back packs printed with Lot No. and expiry date

Product Applications

Fast nonvolumetric transfer of liquids

Drop counting

Blood banks

Hematology

Serology

Bacteriology

Immunology

Chemistry

General laboratory applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	E1414-1311
Bulb draw	3.5 ml
Material	Polyethylene
Pack Size	500 Individually wrapped sachets (1 Pack × 500 Individually wrapped sachets)
Autoclavable	No
Sterile	Yes
Graduations	No
Diameter	4.4 mm
Length	150 mm
Volume	3 ml

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about 3.0 ml Fine Tip Pasteur Pipette (Sterile), Ind. Wrapped

A range of polyethylene pipettes which are manufactured to ISO standards. Ideal for aliquots, pasteur pipettes transfer small amounts of liquid with precision

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de