



Aluminium Sealing Film, 60 µm (PCR)



Product Highlights

- Suitable for standard PCR
- Thicker seal for greater security
- Suitable for short-term storage
- Easy to pierce with a needle, metal probe or pipette tip
- For use with polypropylene, polystyrene, polycarbonate and polyethylene plates
- Suitable for use between -80 °C and 120 °C

Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	E2796-0792
Overall Dimensions	146mm x 79mm
Working Dimensions	119 mm x 79 mm
Sterile	No
Pack Size	100 Pcs. (1 Box x 100 Pcs.)
Material	Aluminium
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
Pyrogen free	No
RNase free	Yes
Color	Silver
Plate type compatibility	Polypropylene, Polystyrene, Polycarbonate, Polyethylene
Max. temp.	120 °C
Temperature range	-80 to 120 °C
Temperature min. (SL)	-80 °C
Breathable	no
Easy to Pierce with Needle or Metal Probe	yes
Easy to Pierce with Pipette Tip	yes
Low 'Tack to Touch' Adhesive	no
Low Auto-Fluorescence	no
Opaque	yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.


Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



Optically clear	no
Resistant to DMSO	no
Seal	Self-adhesive
Short-Term Storage and Incubation	yes
Suitable for PCR	yes
Suitable for Real-Time PCR (qPCR)	no
Superior Optical Clarity	no

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.


Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about Aluminium Sealing Film, 60 µm (PCR)

Aluminium Sealing Film, 60 µm, with acrylate adhesive. The seal is piercable and has high solvent resistance. Can be used for standard PCR, short and long-term storage, including transportation at low temperatures. Use for the storage of organic solvents, acids, and alkalines.

Suitable for use with all plate material types. The use of a seal applicator is advised; strength of sealing is dependant on application pressure.



Accessories

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 Plate Seal Applicator	5 Pcs. (1 Pack × 5 Pcs.)	I2928-7355
 Microtitre Sealing Brayer	1 Piece (1 Pack × 1 Piece)	E9127-2940

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de