



Polyester PCR Sealing Film, Clear



Product Highlights

- This PCR sealing film offers high performance for robotic applications and extreme conditions
- Excellent seal minimizes evaporation and well-to-well cross contamination
- Fit-to-plate dimensions and chamfered corners result in no plate overhang for robotic applications
- Temperature range: - 40°C to + 120°C
- Compatible with thermal cyclers
- Certified RNase, DNase and nucleic acid free

Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	E2796-0100
Overall Dimensions	141 x 78 mm
Working Dimensions	119 x 78 mm
Pack Size	100 Pcs. (1 Box x 100 Pcs.)
Material	Polyester
DNA free	No
DNase free	Yes
Pyrogen free	No
RNase free	Yes
Color	Clear
Plate type compatibility	Polypropylene, Polystyrene, Polycarbonate, Polyethylene
Max. temp.	120 °C
Temperature range	-40 to 120 °C
Temperature min. (SL)	-40 °C
Easy to Pierce with Needle or Metal Probe	no
Easy to Pierce with Pipette Tip	no
Low 'Tack to Touch' Adhesive	no
Low Auto-Fluorescence	no
Opaque	no
Optically clear	yes
Resistant to DMSO	no
Seal	Self-adhesive

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



Short-Term Storage and Incubation	yes
Suitable for PCR	yes
Suitable for Real-Time PCR (qPCR)	no

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about Polyester PCR Sealing Film, Clear



Plates seals are quick to apply and to remove. Eight self-adhesive seals offered by STARLAB are suitable for standard PCR. The Xtra-Clear Advanced Polyolefin Seal is also suitable for Real-Time PCR and the seal's superior optical clarity yields consistent results. All seals are suitable for short-term storage and incubation.

Important note for Real-Time PCR (qPCR)

Please check that both the plate AND the sealing option are suitable for Real-Time PCR applications.



Accessories

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 Plate Seal Applicator	5 Pcs. (1 Pack × 5 Pcs.)	I2928-7355
 Microtitre Sealing Brayer	1 Piece (1 Pack × 1 Piece)	E9127-2940

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de