



T-75 CytoOne® Flask, Non-Treated



Product Highlights

- All flasks 100% pressure tested for leak-free assurance
- Stacking rims with vented design for thermal equalisation
- Anti-tilt skirts and base bars improve stability
- Wide necks provide easy access while minimising contamination
- Gamma sterilised (SAL 10^{-5}) and certified RNase, DNase, DNA and Pyrogen free
- Optical clarity provides distortion-free microscopy
- Easy open, resealable zip top bags protect unused flasks

Product Applications

Expansion and culture of adherent and suspension cells
Performance of cell-based assays and microscopic analysis

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	CC7672-4175
Pack Size	100 Flasks (20 Bags × 5 Flasks)
Growth area	75 cm²
Sterile	Yes
Max. volume with neck	290 ml
Max. volume without neck	279 ml
Working volume	20 ml
Flask Type	T-75
Cap style	Two-position plug seal caps
DNA free	Yes
DNase free	Yes
Pyrogen free	Yes
RNase free	Yes
Color	White
Neck	Canted
Non-cytotoxic	Yes
Treatment	Non-Treated

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de

















More informations about T-75 CytoOne® Flask, Non-Treated

Pressure tested

Choose non-treated for hybridoma and suspension cultures. The two-position plug seal cap lets you choose secure manual venting or an airtight seal. Flasks feature convenient moulded volume graduations and frosted writing surfaces. The T-25, T-75, and T-150 flasks have traceable lot numbers etched into the base.





All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	T-25 CytoOne® Flask, TC-Treated Flask Type: T-25 Treatment: TC Treated	300 Flasks (30 Bags × 10 Flasks)	CC7682-4325
	 Create PDF Data Sheet		
	T-75 CytoOne® Flask, TC-Treated Flask Type: T-75 Treatment: TC Treated	100 Flasks (20 Bags × 5 Flasks)	CC7682-4175
	 Create PDF Data Sheet		
	T-150 CytoOne® Flask, TC-Treated Flask Type: T-150 Treatment: TC Treated	40 Flasks (8 Bags × 5 Flasks)	CC7682-4415
	 Create PDF Data Sheet		
	T-225 CytoOne® Flask, TC-Treated Flask Type: T-225 Treatment: TC Treated	25 Flasks (5 Bags × 5 Flasks)	CC7682-4225
	 Create PDF Data Sheet		
	T-25 CytoOne® Flask, Non-Treated Flask Type: T-25 Treatment: Non-Treated	300 Flasks (30 Bags × 10 Flasks)	CC7672-4325
	 Create PDF Data Sheet		
	T-75 CytoOne® Flask, Non-Treated Flask Type: T-75 Treatment: Non-Treated	100 Flasks (20 Bags × 5 Flasks)	CC7672-4175
	 Create PDF Data Sheet		
	T-150 CytoOne® Flask, Non-Treated Flask Type: T-150 Treatment: Non-Treated	40 Flasks (8 Bags × 5 Flasks)	CC7672-4415
	 Create PDF Data Sheet		

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Hölftigbaum 38
22143 Hamburg
Email: info@starlab.de



PRODUCT NAME		PACKAGING SIZE	ART. NO.
	T-225 CytoOne® Flask, Non-Treated	25 Flasks (5 Bags × 5 Flasks)	CC7672-4225
	Flask Type: T-225 Treatment: Non-Treated		
	 Create PDF Data Sheet		

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de