



T-75 CytoOne® Flask, TC-Treated, Vented



Product Highlights

- All flasks 100% pressure tested for leak-free assurance
- Stacking rims with vented design for thermal equalisation
- Anti-tilt skirts and base bars improve stability
- Wide necks provide easy access while minimising contamination
- Gamma sterilised (SAL 10^{-5}) and certified RNase, DNase, DNA and Pyrogen free
- Optical clarity provides distortion-free microscopy
- Easy open, resealable zip top bags protect unused flasks

Product Applications

Expansion and culture of adherent and suspension cells
Performance of cell-based assays and microscopic analysis

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	CC7682-4875
Growth area	75 cm ²
Sterile	Yes
Max. volume with neck	290 ml
Max. volume without neck	279 ml
Working volume	20 ml
Pack Size	100 Flasks (20 Bags × 5 Flasks)
Flask Type	T-75
Cap style	Vented filter caps
DNA free	Yes
DNase free	Yes
Pyrogen free	Yes
RNase free	Yes
Color	Green
Neck	Canted
Non-cytotoxic	Yes
Treatment	TC Treated

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about T-75 CytoOne® Flask, TC-Treated, Vented















Leak-free assurance

Tissue culture treated for optimum cell attachment and growth or non-treated for hybridoma and suspension cultures. The vented filter caps have a 0.2 µm hydrophobic membrane to exclude bacteria and fungal contaminants. Flasks feature convenient moulded volume graduations and frosted writing surfaces. The T-25, T-75, and T-150 flasks have traceable lot numbers etched into the base.

- › All flasks 100% pressure tested for leak-free assurance
- › Stacking rims with vented design for thermal equalisation
- › Anti-tilt skirts and base bars improve stability
- › Wide necks provide easy access while minimising contamination
- › Optical clarity provides distortion-free microscopy
- › Convenient T150 size for easier sub-culturing and upscaling
- › Easy open, resealable zip-top bags protect unused flasks
- › Gamma sterilised
- › Certified RNase, DNase, DNA and pyrogen free




All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 T-25 CytoOne® Flask, TC-Treated, Vented Flask Type: T-25 Treatment: TC Treated  Create PDF Data Sheet	300 Flasks (30 Bags × 10 Flasks)	CC7682-4825
 T-75 CytoOne® Flask, TC-Treated, Vented Flask Type: T-75 Treatment: TC Treated  Create PDF Data Sheet	100 Flasks (20 Bags × 5 Flasks)	CC7682-4875
 T-150 CytoOne® Flask, TC-Treated, Vented Flask Type: T-150 Treatment: TC Treated  Create PDF Data Sheet	40 Flasks (8 Bags × 5 Flasks)	CC7682-4815
 T-225 CytoOne® Flask, TC-Treated, Vented Flask Type: T-225 Treatment: TC Treated  Create PDF Data Sheet	25 Flasks (5 Bags × 5 Flasks)	CC7682-4822
 T-25 CytoOne® Non-Treated Flasks with Vented Filter Caps Flask Type: T-25 Treatment: Non-Treated  Create PDF Data Sheet	300 Flasks (30 Bags × 10 Flasks)	CC7672-4825
 T-75 CytoOne® Non-Treated Flasks with Vented Filter Caps Flask Type: T-75 Treatment: Non-Treated  Create PDF Data Sheet	100 Flasks (20 Bags × 5 Flasks)	CC7672-4875
 T-150 CytoOne® Non-Treated Flasks with Vented Filter Caps Flask Type: T-150 Treatment: Non-Treated  Create PDF Data Sheet	40 Flasks (8 Bags × 5 Flasks)	CC7672-4815

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <p>T-225 CytoOne® Non-Treated Flasks with Vented Filter Caps Flask Type: T-225 Treatment: Non-Treated</p> <p> Create PDF Data Sheet</p>	25 Flasks (5 Bags × 5 Flasks)	CC7672-4822

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de