



70 µm Cell Strainer, White (Sterile)



Product Highlights

- Uniform nylon mesh size of 70 µm (colour-coded)
- Extended lip on the strainer enables aseptic handling with forceps
- Designed to fit perfectly into a 50 ml centrifuge tube
- RNase, DNase and Pyrogen free
- Gamma sterilised (SAL 10⁻⁶)

Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	CC81111-0072
Pack Size	50 Pcs. (50 Box × 1 Pcs.)
Non-cytotoxic	Yes
Pyrogen free	Yes
RNase free	Yes
DNase free	Yes
Sterile	Yes
Pore size	70 µm
Color	White

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about 70 µm Cell Strainer, White (Sterile)







STARLAB Cell Strainers are intended as a fast, easy and simple alternative to traditional filtration when dissociating stem cells and other tissue-derived primary cells. They yield consistently uniform single-cell suspensions that are ideal for preparation of samples for flow cytometry, Fluorescence-Activated Cell Sorting (FACS), and other applications related to cell separation. Available in 40 µm, 70 µm, and 100 µm color-coded mesh variants. Cell strainers are made from USP Class VI materials.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	40 µm Cell Strainer, Blue (Sterile) Pore size: 40 µm  Create PDF Data Sheet	50 Pcs. (50 Box × 1 Pcs.)	CC8111-0042
	70 µm Cell Strainer, White (Sterile) Pore size: 70 µm  Create PDF Data Sheet	50 Pcs. (50 Box × 1 Pcs.)	CC8111-0072
	100 µm Cell Strainer, Yellow (Sterile) Pore size: 100 µm  Create PDF Data Sheet	50 Pcs. (50 Box × 1 Pcs.)	CC8111-0102

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de