



M210-LAB Printer



Product Highlights

- Label materials for laboratory applications ranging from Liquid Nitrogen (-196 °C) to Autoclave (121 °C)
- 133 symbols: Greek, arrows, warning signs and more
- Barcoding
- Automated serialisation
- Time and date stamp
- Create labels easily with on-board label design wizards
- Save up to 12 label styles for future use
- Auto label function cuts label according to vial size
- Auto shut-off saves battery life
- 2-year warranty
- Print resolution: 203 dpi
- Max. labels per day: 100

Product Applications

Label printing

Plate labelling

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

| | |
|------------------------|--|
| Art. No. | N3000-2000 |
| Pack Size | 1 Piece (1 Box × 1 Piece) |
| Dimensions (W × D × H) | 115 mm × 34 mm × 242 mm |
| Dimensions (W × D) | 115 mm × 34 mm |
| Max. label width [mm] | 19 mm |
| Max. print width [mm] | 914 mm |
| Weight | 0.7 kg |
| Power | 6 x AA Batteries (not included) or AC (optional adapter) |
| Agency approval | Printer: FCC Class B Approved; Labels: UL and CSA Approved |
| Barcoding | Yes (Code 128. Code 39. Data Matrix. QR Code) |
| Display | LCD |
| Keyboard | ABC |
| Memory | Stores 12 files |
| Print resolution [dpi] | 203 |
| Print speed | 10 mm per second rpm |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about M210-LAB Printer

Tough, portable printing

The Brady M210-LAB portable printer has a tough exterior and smart printing capabilities. High quality label materials designed specifically for the laboratory environment create durable, professional looking labels for your samples. An array of easy to load cartridges are available, to cover a multitude of laboratory applications such as liquid nitrogen, freezing, boiling and autoclave as well as general purpose labelling.



Accessories

| PRODUCT NAME | | PACKAGING SIZE | ART. NO. |
|---|--|----------------------------|------------|
|  | BMP21™ AC Adaptor, continental | 1 Piece (1 Pack × 1 Piece) | N3000-2002 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de