



## M210-LAB Printer



### Product Highlights

- Label materials for laboratory applications ranging from Liquid Nitrogen (-196 °C) to Autoclave (121 °C)
- 133 symbols: Greek, arrows, warning signs and more
- Barcoding
- Automated serialisation
- Time and date stamp
- Create labels easily with on-board label design wizards
- Save up to 12 label styles for future use
- Auto label function cuts label according to vial size
- Auto shut-off saves battery life
- 2-year warranty
- Print resolution: 203 dpi
- Max. labels per day: 100

## Product Applications

Label printing

Plate labelling

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	N3000-2000
Pack Size	1 Piece (1 Box × 1 Piece)
Dimensions (W × D × H)	115 mm × 34 mm × 242 mm
Dimensions (W × D)	115 mm × 34 mm
External width	115 mm
Max. label width [mm]	19 mm
Max. print width [mm]	914 mm
Weight	0.7 kg
Power	6 x AA Batteries (not included) or AC (optional adapter)
Agency approval	Printer: FCC Class B Approved; Labels: UL and CSA Approved
Barcoding	Yes (Code 128, Code 39, Data Matrix, QR Code)
Display	LCD
Keyboard	ABC
Memory	Stores 12 files
Print resolution [dpi]	203
Print speed	10 mm per second rpm

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about M210-LAB Printer

### **Tough, portable printing**

The Brady M210-LAB portable printer has a tough exterior and smart printing capabilities. High quality label materials designed specifically for the laboratory environment create durable, professional looking labels for your samples. An array of easy to load cartridges are available, to cover a multitude of laboratory applications such as liquid nitrogen, freezing, boiling and autoclave as well as general purpose labelling.



## Accessories

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	<a href="#">BMP21™ AC Adaptor, continental</a>	1 Piece (1 Pack × 1 Piece)	N3000-2002

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)