



1,250 µl XL graduated TipOne® Tip, Natural, Bag (non-sterile)



Product Highlights

- Graduated at 100 µl, 250 µl, 500 µl and 1000 µl
- Universal tip fit
- High quality virgin polypropylene tips, with extremely smooth surfaces for reduced sample retention
- Certified RNase, DNase, DNA and Pyrogen free.

Product Applications

Pipetting liquids

Dispensing liquids

Mixing liquids

Filling plates and reaction vessels

Phase extraction and removal of supernatants

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	S1112-1020
Sterile	No
Autoclavable	Yes
Filter tip	Non-filter
Volume range	100 - 1,250 µl
Volume	1250 µl
Pack Size	1,000 Tips (1 Bag × 1000 Tips)
Color Code Wafer	Blue
Tip Type	XL
Material	Polypropylene
Graduations	Yes
Color	Natural
For use with brand + product (SL)	<p>Abimed: BioPette-E, Discovery Comfort, LABMATE Optima</p> <p>Biohit: ePet, mLine, Proline</p> <p>BRAND: Transferpette, Transferpette Electronic, Transferpette Pro, Transferpette S</p> <p>Capp: CappAero Comfort, CappAero Universal, Capptronic Single, ecopipette, Singlechannel</p> <p>Eppendorf: Reference, Reference 2, Research, Research 3 neo, Research plus, Research pro, Xplorer</p> <p>Gilson: Pipetman Concept, Pipetman L, Pipetman Neo, Pipetman P, Pipetman Ultra</p> <p>Hamilton: SoftGrip Pipette, SofTouch</p> <p>HTL: Clinipet+, Discovery Comfort, Discovery Pro, Labmate Pro, Optipette</p> <p>Nichiryo: 31N, 50F, 50V, NLE, NPLO, NPX</p> <p>Rainin: Classic, EDP1, EDP-3 Plus, Pipet-Lite, Pipet-Plus</p> <p>Sartorius: ePet, mLine, Proline, Proline Plus, Tacta</p> <p>Socorex: Acura 815 manual, Acura 825, Acura electro 925, Calibra digital</p> <p>STARLAB: ErgoOne, ErgoOne E, MicroOne, StarPet, StarPet-E</p> <p>Thermo Fisher Scientific Inc.: Finn Focus, Finnpipette Digital, Finnpipette F1/F2, Finnpipette Novus</p>

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
 Neuer Höltigbaum 38
 22143 Hamburg
 Email: info@starlab.de















More informations about 1,250 µl XL graduated TipOne® Tip, Natural, Bag (non-sterile)

High quality, reliable [pipette tips](#) with a universal tip fit – [TipOne](#)® has it all! If you are looking for accuracy and precision you can trust, [TipOne](#)® is the one tip for you.

Loose tips are supplied in a convenient, resealable bag. The tamper-proof perforation strip on the bag provides easy opening, whilst the resealable zip-lock provides contamination-free storage of opened product.











All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <p>1,250 µl XL graduated TipOne® Tip, Natural, Racks (non-sterile) Sterile: No Pack Size: 960 Tips (10 Racks × 96 Tips)</p> <p> Create PDF Data Sheet</p>	960 Tips (10 Racks × 96 Tips)	S1112-1820
 <p>1,250 µl XL graduated TipOne® Tip, Natural, Racks (non-sterile), Case Sterile: No Pack Size: 3,840 Tips (40 Racks × 96 Tips)</p> <p> Create PDF Data Sheet</p>	3,840 Tips (40 Racks × 96 Tips)	S1112-1820-C
 <p>1,250 µl XL graduated TipOne® Tip, Natural, Racks (sterile) Sterile: Yes Pack Size: 960 Tips (10 Racks × 96 Tips)</p> <p> Create PDF Data Sheet</p>	960 Tips (10 Racks × 96 Tips)	S1112-1830
 <p>1,250 µl XL graduated TipOne® Tip, Natural, Racks (sterile), Case Sterile: Yes Pack Size: 3,840 Tips (40 Racks × 96 Tips)</p> <p> Create PDF Data Sheet</p>	3,840 Tips (40 Racks × 96 Tips)	S1112-1830-C
 <p>1,250 µl XL graduated TipOne® Tip, Natural, Refills (non-sterile) Sterile: No Pack Size: 960 Tips (10 Refills × 96 Tips)</p> <p> Create PDF Data Sheet</p>	960 Tips (10 Refills × 96 Tips)	S1112-1720
 <p>1,250 µl XL graduated TipOne® Tip, Natural, Refills (non-sterile), Case Sterile: No Pack Size: 9,600 Tips (100 Refills × 96 Tips)</p> <p> Create PDF Data Sheet</p>	9,600 Tips (100 Refills × 96 Tips)	S1112-1720-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
 Neuer Höltigbaum 38
 22143 Hamburg
 Email: info@starlab.de



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 1,250 µl XL graduated TipOne® Tip, Natural, Refills (sterile) Sterile: Yes Pack Size: 960 Tips (10 Refills × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Refills × 96 Tips)	S1112-1730
 1,250 µl XL graduated TipOne® Tip, Natural, Refills (sterile), Case Sterile: Yes Pack Size: 3,840 Tips (40 Refills × 96 Tips)  Create PDF Data Sheet	3,840 Tips (40 Refills × 96 Tips)	S1112-1730-C
 1,250 µl XL graduated TipOne® Tip, Natural, Bag (non-sterile) Sterile: No Pack Size: 1,000 Tips (1 Bag × 1000 Tips)  Create PDF Data Sheet	1,000 Tips (1 Bag × 1000 Tips)	S1112-1020
 1,250 µl XL graduated TipOne® Tip, Natural, Bags (non-sterile), Case Sterile: No Pack Size: 16,000 Tips (16 Bags × 1000 Tips)  Create PDF Data Sheet	16,000 Tips (16 Bags × 1000 Tips)	S1112-1020-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de