



300 µl graduated TipOne® Tip, Natural, Bag (non-sterile)



Product Highlights

- Graduated at 100 µl, 200 µl and 300 µl
- Universal fit. One tip fits all! The tip collar is engineered to fit all popular pipettes.
- Highly polished moulds for an extremely smooth surface to reduce sample retention
- Highest quality product. Tips manufactured from pure, virgin polypropylene to ensure the highest quality product.
- Certified RNase, DNase, DNA and Pyrogen free.

Product Applications

Pipetting liquids

Dispensing liquids

Mixing liquids

Filling plates and reaction vessels

Phase extraction and removal of supernatants

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	S1110-9000
Pack Size	1,000 Tips (1 Bag × 1000 Tips)
Color Code Wafer	Green
Sterile	No
Autoclavable	Yes
Filter tip	Non-filter
Volume range	1.0 - 300 µl
Volume	300 µl
Tip Type	standard
Material	Polypropylene
Graduations	Yes
Color	Natural
For use with	<p>Abimed: BioPette-E, Discovery Comfort, LABMATE Optima</p> <p>Biohit®: eLine, ePet, mLine, Picus, Proline</p> <p>BRAND®: Transferpette, Transferpette Electronic, Transferpette Pro, Transferpette S</p> <p>Capp®: CappAero Comfort, CappAero Multi, CappAero Universal, Capptronic Multi, Capptronic Single, ecopipette, Singlechannel</p> <p>Eppendorf®: Reference, Reference 2, Research, Research plus, Research pro, Xplorer</p> <p>Gilson®: Pipetman Concept, Pipetman L, Pipetman Neo, Pipetman P, Pipetman Ultra</p> <p>Hamilton®: SoftGrip Pipette, SofTouch</p> <p>HTL: Clinipet+, Discovery Comfort, Discovery Pro, Labmate Pro, Optipette</p> <p>Nichiryo®: 00-NP7, 50F, 50V, NLE, NPLO, NPX</p> <p>Rainin®: Classic, EDP1, EDP-3 Plus, Pipet-Lite, Pipet-Plus</p> <p>Sartorius®: eLine, ePet, mLine, Picus, Picus NxT, Proline, Proline Plus, Tacta</p> <p>Socorex®: Acura 815 manual, Acura 825, Acura 855, Acura electro 925, Acura electro 955, Calibra digital</p> <p>STARLAB®: ErgoOne, ErgoOne E, MicroOne, StarPet, StarPet-E</p> <p>Thermo Fisher Scientific Inc.: Finn Focus, Finn Focus Micro, Finn Focus Univ., Finnpipette Digital, Finnpipette F1/F2, Finnpipette Novus, Finnpipette Novus Turquoise, Finnpipette Novus Yellow</p>

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about 300 µl graduated TipOne® Tip, Natural, Bag (non-sterile)













Loose tips are supplied in a convenient, self-standing bag which provides easy access and storage. The tamper-proof perforation strip on the bag enables easy opening. A resealable zip-lock provides contamination-free storage of opened product. The bag is made from recyclable polyethylene.

UltraPoint orifice

Thin-wall orifice helps prevent sample hang up and makes it easier to touch off the last drop of sample











All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 300 µl graduated TipOne® Tip, Natural, Racks (non-sterile) Sterile: No Pack Size: 960 Tips (10 Racks × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Racks × 96 Tips)	S1110-9800
 300 µl graduated TipOne® Tip, Natural, Racks (non-sterile), Case Sterile: No Pack Size: 7,680 Tips (80 Racks × 96 Tips)  Create PDF Data Sheet	7,680 Tips (80 Racks × 96 Tips)	S1110-9800-C
 300 µl graduated TipOne® Tip, Natural, Racks (sterile) Sterile: Yes Pack Size: 960 Tips (10 Racks × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Racks × 96 Tips)	S1110-9810
 300 µl graduated TipOne® Tip, Natural, Racks (sterile), Case Sterile: Yes Pack Size: 7,680 Tips (80 Racks × 96 Tips)  Create PDF Data Sheet	7,680 Tips (80 Racks × 96 Tips)	S1110-9810-C
 300 µl graduated TipOne® Tip, Natural, Refills (non-sterile) Sterile: No Pack Size: 960 Tips (10 Refills × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Refills × 96 Tips)	S1110-9700
 300 µl graduated TipOne® Tip, Natural, Refills (non-sterile), Case Sterile: No Pack Size: 9,600 Tips (100 Refills × 96 Tips)  Create PDF Data Sheet	9,600 Tips (100 Refills × 96 Tips)	S1110-9700-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 300 µl graduated TipOne® Tip, Refills (sterile) Sterile: Yes Pack Size: 960 Tips (10 Refills × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Refills × 96 Tips)	S1110-9710
 300 µl graduated TipOne® Tip, Natural, Refills (sterile), Case Sterile: Yes Pack Size: 7,680 Tips (80 Refills × 96 Tips)  Create PDF Data Sheet	7,680 Tips (80 Refills × 96 Tips)	S1110-9710-C
 300 µl graduated TipOne® Tip, Natural, Bag (non-sterile) Sterile: No Pack Size: 1,000 Tips (1 Bag × 1000 Tips)  Create PDF Data Sheet	1,000 Tips (1 Bag × 1000 Tips)	S1110-9000
 300 µl graduated TipOne® Tip, Natural, Bags (non-sterile), Case Sterile: No Pack Size: 25,000 Tips (25 Bags × 1000 Tips)  Create PDF Data Sheet	25,000 Tips (25 Bags × 1000 Tips)	S1110-9000-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de