



## ID Tape Dispenser



### Product Highlights

- Custom-made acrylic dispenser which holds both lengths of tape
- Suction feet ensure the dispenser stays firmly on the bench, even when using one hand to dispense the tape
- Core inserts for 55 m long tape available separately

## Product Applications

Labelling

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

|                        |                           |
|------------------------|---------------------------|
| Art. No.               | E9026-0000                |
| Dimensions (W × D × H) | 205 mm × 185 mm × 125 mm  |
| Dimensions (W × D)     | 205 mm × 185 mm           |
| Pack Size              | 1 Piece (1 Box × 1 Piece) |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about ID Tape Dispenser


### Dispense with disorganisation!

- > Custom-made acrylic [dispenser](#) which holds both lengths of tape
- > Suction feet ensures the [dispenser](#) stays firmly on the bench, even when using one hand to dispense the tape
- > Core inserts to hold the 55 m long tape in the [dispenser](#) available separately

[Dispenser](#) dimensions (WxDxH): 205 x 185 x 125 mm



## Accessories

| PRODUCT NAME  | PACKAGING SIZE          | ART. NO.   |
|---|-------------------------|------------|
|  <a href="#">Starlab ID Tape Core Inserts, Pack of 5</a> | 5 Pcs. (1 Box x 5 Pcs.) | E9026-2005 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)