



## 1,000 µl XL graduated TipOne® Filter Tip, Natural, Racks (sterile)



### Product Highlights

- Graduated at 100 µl, 250 µl, 500 µl and 1000 µl
- Highly effective HDPE filter that contains no additives to interfere with samples
- Universal tip fit
- High-quality virgin polypropylene tips, with extremely smooth surfaces for reduced sample retention
- Sterile and Certified RNase, DNase, DNA and Pyrogen free
- Racks can be retained and reused with TipOne® refills
- All components, excluding tips and shrink wrap, are recyclable

## Product Applications

DNA applications (e.g. PCR)

RNA applications (e.g. Gene expression analysis)

Protein applications (e.g. Antibody Research)

Cell Culture applications (e.g. Media)

Applications with radio actives

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	S1122-1830
Pack Size	960 Tips (10 Racks × 96 Tips)
Color Code Wafer	Blue
Sterile	Yes
Autoclavable	No
Filter tip	Filter
Volume range	100 - 1,000 µl
Volume	1000 µl
Tip Type	XL
Material	Polypropylene
Graduations	Yes
Color	Natural
Sustainability attributes	Racks can be retained and reused with TipOne® refills. Rack base made with ~25 % recycled plastic. All components, excluding tips and shrink wrap, are recyclable.
For use with	<b>Abimed</b> : BioPette-E, Discovery Comfort, LABMATE Optima <b>Biohit</b> ®: ePet, mLine, Proline <b>BRAND</b> ®: Transferpette, Transferpette Electronic, Transferpette Pro, Transferpette S <b>Capp</b> ®: CappAero Comfort, CappAero Universal, Capptronic Single, ecopipette, Singlechannel <b>Eppendorf</b> ®: Reference, Reference 2, Research, Research plus, Research pro, Xplorer <b>Gilson</b> ®: Pipetman L, Pipetman Neo, Pipetman P, Pipetman Ultra <b>Hamilton</b> ®: SoftGrip Pipette, SofTouch <b>HTL</b> : Clinipet+, Discovery Comfort, Discovery Pro, Labmate Pro, Optipette <b>Nichiryo</b> ®: 31N, 50F, 50V, NLE, NPLO, NPX <b>Rainin</b> ®: Classic, EDP1, EDP-3 Plus, Pipet-Lite, Pipet-Plus <b>Sartorius</b> ®: ePet, mLine, Proline, Proline Plus, Tacta <b>Socorex</b> ®: Acura 815 manual, Acura 825, Acura electro 925, Calibra digital <b>STARLAB</b> ®: ErgoOne, ErgoOne E, MicroOne, StarPet, StarPet-E <b>Thermo Fisher Scientific Inc.</b> : Finn Focus, Finnpipette Digital, Finnpipette F1/F2, Finnpipette Novus

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about 1,000 µl XL graduated TipOne® Filter Tip, Natural, Racks (sterile)

### Nothing goes through!

Protect your samples from unwanted contamination with TipOne® filter tips. Our high-quality filters prevent aerosol cross-contamination without trapping your valuable samples. Plus, they contain no additives that could interfere with your samples.

You can't always see fine aerosols entering your pipette, but they can contaminate other samples later. TipOne® filter tips block over 99% of aerosols\*, making them one of the best filters available.

Sterile TipOne® tips are validated according to EN ISO 11137-2:2007.

### Join our journey to sustainable pipetting.









TipOne® tip racks are built to last. When you have used your filter tips, the rack can be reused again and again with the cost-effective TipOne® refills (non-sterile, sterile and sterile filter tip options available). TipOne® racks have been tested and they can withstand at least 100 autoclave cycles (rack with wafer inside) without damage.

The base of the rack is made with 25 % recycled plastic. All components, excluding tips and shrink, are recyclable. Sign up for our TipOne® recycling service!

\*See TipOne® brochure in downloads section for details.



## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">1,000 µl XL graduated TipOne® Filter Tip, Natural, Refills (sterile)</a> <b>Pack Size:</b> 960 Tips (10 Refills × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (10 Refills × 96 Tips)	S1122-1730
 <a href="#">1,000 µl XL graduated TipOne® Filter Tip, Natural, Refills (sterile) Case</a> <b>Pack Size:</b> 3,840 Tips (40 Refills × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	3,840 Tips (40 Refills × 96 Tips)	S1122-1730-C
 <a href="#">1,000 µl XL graduated TipOne® Filter Tip, Natural, Racks (sterile)</a> <b>Pack Size:</b> 960 Tips (10 Racks × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (10 Racks × 96 Tips)	S1122-1830
 <a href="#">1,000 µl XL graduated TipOne® Filter Tip, Natural, Racks (sterile) Case</a> <b>Pack Size:</b> 3,840 Tips (40 Racks × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	3,840 Tips (40 Racks × 96 Tips)	S1122-1830-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)