




## Extended Length Filter Tips



### Product Highlights

- Compatible with a wide range of pipettes
- Made from polypropylene
- For some tips sterile racks or bags available

## Product Applications

Pipetting liquids from and into long tubes  
Phase extraction and removal of lower phases

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	See variations
Filter tip	Filter
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Pyrogen free	Yes
RNase free	Yes
Color	Natural

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about Extended Length Filter Tips











Extra long tips for easy handling of liquids in long tubes. Phase extraction and removal of lower phases.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">1.25 ml XL Graduated Filter Tip, Racked (sterile)</a> <b>Volume:</b> 1.25 ml <b>Sterile:</b> Yes <b>Pack Size:</b> 768 Tips (8 Racks × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	768 Tips (8 Racks × 96 Tips)	I1011-3010
 <a href="#">200 µl XL Filter Tip, sterile, Racked</a> <b>Volume:</b> 200 µl <b>Sterile:</b> Yes <b>Pack Size:</b> 800 Tips (8 Racks × 100 Tips)  <a href="#">Create PDF Data Sheet</a>	800 Tips (8 Racks × 100 Tips)	I1011-5830
 <a href="#">200 µl XL Filter Tip, sterile, Racked</a> <b>Volume:</b> 200 µl <b>Sterile:</b> Yes <b>Pack Size:</b> 1,632 Tips (8 Racks × 204 Tips)  <a href="#">Create PDF Data Sheet</a>	1,632 Tips (8 Racks × 204 Tips)	I1011-0830
 <a href="#">1,000 µl XL Bevelled Filter Tip, sterile, Racked</a> <b>Volume:</b> 1000 µl <b>Sterile:</b> Yes <b>Pack Size:</b> 1,000 Tips (10 Racks × 100 Tips)  <a href="#">Create PDF Data Sheet</a>	1,000 Tips (10 Racks × 100 Tips)	I1026-7810
 <a href="#">1,000 µl XL Bevelled Filter Tip, Loose / 1000 tips (1 bag x 1000 tips)</a> <b>Volume:</b> 1000 µl <b>Sterile:</b> No <b>Pack Size:</b> 1,000 Tips (1 Bag × 1000 Tips)  <a href="#">Create PDF Data Sheet</a>	1,000 Tips (1 Bag × 1000 Tips)	E1011-1008

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">200 µl XL Filter Tip, sterile, Ind-Wrap/ 100 tips (100 x individually wrapped sachets)</a> <b>Volume:</b> 200 µl <b>Sterile:</b> Yes <b>Pack Size:</b> 100 Individually wrapped sachets (1 Bag × 100 Individually wrapped sachets)  <a href="#">Create PDF Data Sheet</a>	100 Individually wrapped sachets (1 Bag × 100 Individually wrapped sachets)	I1011-0070
 <a href="#">1,000 µl XL Bevelled Filter Tip</a> <b>Volume:</b> 1000 µl <b>Sterile:</b> Yes <b>Pack Size:</b> 100 Individually wrapped sachets (1 Bag × 100 Individually wrapped sachets)  <a href="#">Create PDF Data Sheet</a>	100 Individually wrapped sachets (1 Bag × 100 Individually wrapped sachets)	I1026-1070

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)