




Flat Gel-Loading Tips



Product Highlights

- Flexible, long and narrow
- Loading samples onto gels (eg. polyacrylamide or agarose)
- Handling of small volumes
- Compatible with a wide range of pipettes
- Different pack types available
- Certified

Product Applications

Gel electrophoresis

Loading samples onto gels (eg. polyacrylamide or agarose)

Handling of small volumes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

| | |
|-----------|------------------------------|
| Art. No. | See variations |
| Sterile | No |
| Volume | 200 µl |
| Pack Size | 200 Tips (1 Rack × 200 Tips) |
| Tip Type | flat gel-loading |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about Flat Gel-Loading Tips





All tips are certified RNase and DNase free. Packaging of these tips differs from the [TipOne](#)[®] range.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



All Variations

| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|--|------------------------------|------------|
|  200 µl Flat Gel-Loading Tip (Tip Thickness 0.20 mm) Tip Thickness: 0.20 mm  Create PDF Data Sheet | 200 Tips (1 Rack × 200 Tips) | I1022-2600 |
|  200 µl Flat Gel-Loading Tip (Tip Thickness 0.40 mm) Tip Thickness: 0.40 mm  Create PDF Data Sheet | 200 Tips (1 Rack × 200 Tips) | I1022-4600 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de