



ErgoOne® Multi-Channel Pipettes



Product Highlights

- Minimal operating forces
- Lightweight ergonomic design
- Colour coded
- Volume lock
- Universal tip fit
- Fully autoclavable

Product Applications

Pipetting liquids

Mixing liquids

Reverse pipetting of "tricky" liquids (e. g. with high viscosity)

Filling of plates, gels, and reaction vessels

Phase extraction and removal of supernatants

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	See variations
Autoclavable	Yes
Operating mode	Manual
Pipetting type	Air-cushion
Operation	Mechanical, separate control button and ejector

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about ErgoOne® Multi-Channel Pipettes

A great tool to speed up your pipetting process

ErgoOne® multi-channel [pipettes](#) are designed in accordance with the latest ergonomic principles. The pipette is extremely easy to use and has minimal operating forces to help prevent hand fatigue. ErgoOne® multi-channel [pipettes](#) provide accurate and precise pipetting in a variety of volume ranges.













As one of the most frequently used instruments in the laboratory, [pipettes](#) must also be comfortable to use and provide the features you need. ErgoOne® [pipettes](#) combine researcher-requested features with the latest ergonomic principles to provide optimum performance. ErgoOne® gives you accuracy and precision you can trust, with a stylish new look you will love.

Lightweight and extremely easy to use, ErgoOne® can be trusted to provide accurate and precise pipetting all day, every day. Use with STARLAB's [TipOne® pipette tips](#) for the perfect pipetting partnership.

ErgoOne® is designed to fit all popular [pipettes](#) without any loss of accuracy and precision.



All Variations





PRODUCT NAME	PACKAGING SIZE	ART. NO.
 0.5 – 10 µl ErgoOne® 8-Channel Pipette Volume range: 0.5 - 10 µl Channel type: 8-channel  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	S7108-0510
 10 – 100 µl ErgoOne® 8-Channel Pipette Volume range: 10 - 100 µl Channel type: 8-channel  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	S7108-1100
 30 – 300 µl ErgoOne® 8-Channel Pipette Volume range: 30 - 300 µl Channel type: 8-channel  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	S7108-3300
 0.5 – 10µl ErgoOne® 12-Channel Pipette Volume range: 0.5 - 10 µl Channel type: 12-channel  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	S7112-0510
 10 – 100µl ErgoOne® 12-Channel Pipette Volume range: 10 - 100 µl Channel type: 12-channel  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	S7112-1100
 30 – 300µl ErgoOne® 12-Channel Pipette Volume range: 30 - 300 µl Channel type: 12-channel  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	S7112-3300

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



Accessories

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 Acrylic Stand for up to 5 Pipettes (Universal)	1 Piece (1 Pack × 1 Piece)	E3099-0000
 Benchtop Support for ErgoOne® Pipettes	5 Pcs. (1 Pack × 5 Pcs.)	E3099-0002
 Self-Adhesive Shelf Holder for one ErgoOne®	1 Piece (1 Pack × 1 Piece)	E3099-0001
 ErgoOne® Carousel Stand for up to 8 Pipettes (Universal)	1 Piece (1 Box × 1 Piece)	S7200-0000

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de