



## ErgoOne® Single-Channel Starter Packs



### Product Highlights

- Three or four ErgoOne® single-channel pipettes
- Three TipOne® racks with compatible tips
- Three or four Self-Adhesive Shelf Holder for one Pipette
- One STARLAB Lab Timer
- Please select the volume range for more detailed info about the content of each Starter Pack.

## Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	See variations
Volume range	2 - 1,000 µl
Volume	1000 µl
Channel type	1-channel
Operating mode	Manual

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about ErgoOne® Single-Channel Starter Packs







Every day, [pipettes](#) are confronted with multiple tasks that must be exact and precise. Reliability is key in the laboratory and your [pipette](#) must be trusted to perform well, day in, day out. Start your new pipetting experience with [ErgoOne](#)® single-channel [pipettes](#) and original [TipOne](#)® tips - the perfect pipetting partnership! Starter Packs also include one Self-Adhesive Shelf Holder for each [pipette](#), plus a STARLAB Lab [Timer](#).

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">ErgoOne® Starter Pack 4 (2.5)</a> Volume range: 0.1 - 1,000 µl  <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	STARTERPACK42.5
 <a href="#">ErgoOne® Starter Pack 4 (10)</a> Volume range: 0.5 - 1,000 µl  <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	STARTERPACK410
 <a href="#">ErgoOne® Starter Pack</a> Volume range: 20 - 1,000 µl  <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	STARTERPACK

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)