



## ErgoOne® E Single-Channel Pipettes



### Product Highlights

- Easy to use, intuitive menu structure
- Lightweight, ergonomic design
- Adjustable finger rest for individual comfort
- Significantly reduced tip loading and ejection forces
- Universal tip fit
- Reliable
- Autoclavable pipette shaft
- 1-year warranty

## Product Applications

Pipetting

Reverse Pipetting of "tricky" liquids (e. g. with high viscosity)

Mixing

Dispense

Electrophoresis (excl. 1000 µl and 5000 µl models)

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	See variations
Autoclavable	Autoclavable pipette shaft
Operating mode	Electronic
Pipetting type	Air-cushion
Operation	Battery operated piston stroke / manual tip ejection

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about ErgoOne® E Single-Channel Pipettes

### Reliable and ergonomic!

ErgoOne® E electronic pipettes provide reliable pipetting every day. Ergonomic and lightweight, ErgoOne® E pipettes are perfect for comfortable pipetting during long periods of use. ErgoOne® E and Starlab's TipOne® pipette tips provide the perfect pipetting partnership, but the optimized design of the tip cone provides the user with a convenient universal tip fit.

Four operating modes provide standard and reverse pipetting, as well as mixing and dispensing. A unique "Electrophoresis" mode (available on all models excluding 1000 µl, 1250 µl, and 5000 µl,) is also available to help you load gels easily. The long-lasting battery can also be charged during use so you are never without your pipette.




## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">ErgoOne® E Single-Channel Pipette Electronic, 0.5 – 10µl</a> <b>Volume range:</b> 0.5 - 10 µl  <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	G9001-0010
 <a href="#">ErgoOne® E Single-Channel Pipette Electronic, 2 – 20µl</a> <b>Volume range:</b> 2.0 - 20 µl  <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	G9001-0020
 <a href="#">ErgoOne® E Single-Channel Pipette Electronic, 10 – 200µl</a> <b>Volume range:</b> 10 - 200 µl  <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	G9001-0200
 <a href="#">ErgoOne® E Single-Channel Pipette Electronic, 50 – 1000µl</a> <b>Volume range:</b> 50 - 1,000 µl  <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	G9001-1000
 <a href="#">ErgoOne® E Single Channel Pipette Electronic, 250 – 5000µl</a> <b>Volume range:</b> 250 - 5,000 µl  <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	G9001-5000

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## Accessories

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">Silicone Oil for Piston (for 5000 µl Single-Channel and all Multi-Channel Pipettes)</a>	1 Piece (1 Pack × 1 Piece)	G9999-1007
 <a href="#">Fluorstatic Grease for Guiding Axis (for Multi-Channel Pipettes)</a>	1 Piece (1 Pack × 1 Piece)	G9999-1008
 <a href="#">Individual Stand (for Single-Channel Pipettes, 10 µl to 1000 µl)</a>	1 Piece (1 Box × 1 Piece)	G9999-1002
 <a href="#">Individual Stand (for 5000 µl Single-Channel Pipette)</a>	1 Piece (1 Box × 1 Piece)	G9999-1003
 <a href="#">Shelf Holder for one ErgoOne® E pipette</a>	1 Piece (1 Box × 1 Piece)	G9999-1009
 <a href="#">Wall Mount for one ErgoOne® E pipette</a>	1 Piece (1 Box × 1 Piece)	G9999-1010
 <a href="#">Universal stand for one ErgoOne® E pipette - all volumes</a>	1 Piece (1 Box × 1 Piece)	G9999-1011
 <a href="#">Replacement Battery for ErgoOne® E</a>	1 Piece (1 Box × 1 Piece)	G9999-1005
 <a href="#">AC Adapter for Europe/UK (230V/50 Hz)</a>	1 Piece (1 Box × 1 Piece)	G9999-1000

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)