



ErgoOne® E Multi-Channel Pipettes



Product Highlights

- Easy to use, intuitive menu structure
- Lightweight, ergonomic design
- Adjustable finger rest for individual comfort
- Significantly reduced tip loading and ejection forces
- Universal tip fit
- Reliable
- Autoclavable pipette shaft
- 12-month warranty

Product Applications

Pipetting

Reverse Pipetting of "tricky" liquids (e. g. with high viscosity)

Mixing

Dispense

Electrophoresis (excl. 1250 µl models)

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

| | |
|----------------|--|
| Art. No. | See variations |
| Autoclavable | Autoclavable pipette shaft |
| Operating mode | Electronic |
| Pipetting type | Air-cushion |
| Operation | Battery operated piston stroke / manual tip ejection |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about ErgoOne® E Multi-Channel Pipettes













Reliable and ergonomic!

ErgoOne® E electronic pipettes provide reliable pipetting every day. Ergonomic and lightweight, ErgoOne® E pipettes are perfect for comfortable pipetting during long periods of use. ErgoOne® E and Starlab's TipOne® pipette tips provide the perfect pipetting partnership, but the optimized design of the tip cone provides the user with a convenient universal tip fit.

Four operating modes provide standard and reverse pipetting, as well as mixing and dispensing. A unique "Electrophoresis" mode (available on all models excluding 1000 µl, 1250 µl, and 5000 µl,) is also available to help you load gels easily. The long-lasting battery can also be charged during use so you are never without your pipette.











All Variations

| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|---|---------------------------|------------|
|  ErgoOne® E 8-Channel Electronic Pipette, 0.5 – 10µl Channel type: 8-channel Volume range: 0.5 - 10 µl  Create PDF Data Sheet | 1 Piece (1 Box × 1 Piece) | G9008-0010 |
|  ErgoOne® E 8-Channel Electronic Pipette, 1.0 – 20µl Channel type: 8-channel Volume range: 1.0 - 20 µl  Create PDF Data Sheet | 1 Piece (1 Box × 1 Piece) | G9008-0020 |
|  ErgoOne® E 8-Channel Electronic Pipette, 5.0 – 100µl Channel type: 8-channel Volume range: 5.0 - 100 µl  Create PDF Data Sheet | 1 Piece (1 Box × 1 Piece) | G9008-0100 |
|  ErgoOne® E 8-Channel Electronic Pipette, 15 – 300µl Channel type: 8-channel Volume range: 15 - 300 µl  Create PDF Data Sheet | 1 Piece (1 Box × 1 Piece) | G9008-0300 |
|  ErgoOne® E 8-Channel Electronic Pipette, 50 – 1250µl Channel type: 8-channel Volume range: 50 - 1,250 µl  Create PDF Data Sheet | 1 Piece (1 Box × 1 Piece) | G9008-1250 |
|  ErgoOne® E 12-Channel Pipette Electronic, 5.0 – 100µl Channel type: 12-channel Volume range: 5.0 - 100 µl  Create PDF Data Sheet | 1 Piece (1 Box × 1 Piece) | G9012-0100 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de











| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|--|---------------------------|------------|
|  ErgoOne® E 12-Channel Electronic Pipette, 0.5 – 10µl Channel type: 12-channel Volume range: 0.5 - 10 µl  Create PDF Data Sheet | 1 Piece (1 Box × 1 Piece) | G9012-0010 |
|  ErgoOne® E 12-Channel Electronic Pipette, 1.0 – 20µl Channel type: 12-channel Volume range: 1.0 - 20 µl  Create PDF Data Sheet | 1 Piece (1 Box × 1 Piece) | G9012-0020 |
|  ErgoOne® E 12-Channel Electronic Pipette, 15 – 300µl Channel type: 12-channel Volume range: 15 - 300 µl  Create PDF Data Sheet | 1 Piece (1 Box × 1 Piece) | G9012-0300 |
|  ErgoOne® E 12-Channel Electronic Pipette, 50 – 1250µl Channel type: 12-channel Volume range: 50 - 1,250 µl  Create PDF Data Sheet | 1 Piece (1 Box × 1 Piece) | G9012-1250 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



Accessories

| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|---|----------------------------|------------|
|  Replacement Battery for ErgoOne® E | 1 Piece (1 Box × 1 Piece) | G9999-1005 |
|  Silicone Oil for Piston (for 5000 µl Single-Channel and all Multi-Channel Pipettes) | 1 Piece (1 Pack × 1 Piece) | G9999-1007 |
|  Fluorstatic Grease for Guiding Axis (for Multi-Channel Pipettes) | 1 Piece (1 Pack × 1 Piece) | G9999-1008 |
|  Individual Stand (for all ErgoOne® E Multi-Channel Pipettes up to 300 µl) | 1 Piece (1 Box × 1 Piece) | G9999-1004 |
|  Shelf Holder for one ErgoOne® E pipette | 1 Piece (1 Box × 1 Piece) | G9999-1009 |
|  Wall Mount for one ErgoOne® E pipette | 1 Piece (1 Box × 1 Piece) | G9999-1010 |
|  Universal stand for one ErgoOne® E pipette - all volumes | 1 Piece (1 Box × 1 Piece) | G9999-1011 |
|  AC Adapter for Europe/UK (230V/50 Hz) | 1 Piece (1 Box × 1 Piece) | G9999-1000 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de