



## StarStore 81 Storage Box



### Product Highlights

- 81-place polypropylene storage box with fixed 9x9 dividers for tubes up to 13mm
- Transparent lift-off lid with frosted writing panel on side
- Vent holes in base to aid draining and lid removal
- Printed reference grid on base
- Temperature use -90 °C to 121 °C
- Dimensions (WxDxH): 130 x 130 x 45 mm
- Max. tube height: 60 mm
- Autoclavable

## Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	See variations
Pack Size	5 Pcs. (1 Pack × 5 Pcs.)
Dimensions (W × D × H)	130 mm × 130 mm × 45 mm
Dimensions (W × D)	130 mm × 130 mm
Max. tube length	60 mm
Autoclavable	Yes
Material	Polypropylene
Color	Various
Places	81

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)




## More informations about StarStore 81 Storage Box

















81-place polypropylene storage box with fixed, 9x9 dividers for tubes up to 13mm, eg. 1.5/2.0 ml microcentrifuge tubes or screw cap microtubes. The non-hinged, lift-off lid enables tubes of varying heights to be stored inside. Standard 130x130mm footprint for use in freezer racking.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	<a href="#">StarStore 81, Natural</a> Color: Natural	5 Pcs. (1 Pack × 5 Pcs.)	I2381-5040
	 <a href="#">Create PDF Data Sheet</a>		
	<a href="#">StarStore 81, Neon Blue</a> Color: Neon blue	5 Pcs. (1 Pack × 5 Pcs.)	I2381-5041
	 <a href="#">Create PDF Data Sheet</a>		
	<a href="#">StarStore 81, Neon Green</a> Color: Neon green	5 Pcs. (1 Pack × 5 Pcs.)	I2381-5042
	 <a href="#">Create PDF Data Sheet</a>		
	<a href="#">StarStore 81, Neon Orange</a> Color: Neon orange	5 Pcs. (1 Pack × 5 Pcs.)	I2381-5043
	 <a href="#">Create PDF Data Sheet</a>		
	<a href="#">StarStore 81, Neon Pink</a> Color: Neon pink	5 Pcs. (1 Pack × 5 Pcs.)	I2381-5044
	 <a href="#">Create PDF Data Sheet</a>		
	<a href="#">StarStore 81, Neon Yellow</a> Color: Neon yellow	5 Pcs. (1 Pack × 5 Pcs.)	I2381-5046
	 <a href="#">Create PDF Data Sheet</a>		
	<a href="#">StarStore 81, Neon Mixed</a> Color: Neon mixed	5 Pcs. (1 Pack × 5 Pcs.)	I2381-5048
	 <a href="#">Create PDF Data Sheet</a>		
	<a href="#">StarStore 81, Black</a> Color: Black	5 Pcs. (1 Pack × 5 Pcs.)	I2381-5049
	 <a href="#">Create PDF Data Sheet</a>		

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)