



StarStore 100 Storage Box



Product Highlights

- 100-place polypropylene storage box with fixed 10x10 dividers for tubes up to 13mm
- Hinged, lockable lid with positioning lugs for stable stacking
- Moulded reference grid on base.
- Temperature use -90°C to 121°C
- Autoclavable
- Dimensions (WxDxH): 140 x 155 x 55 mm

Product Applications

Sample storage

Storage of microcentrifuge tubes or screw cap microtubes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	See variations
Pack Size	5 Pcs. (1 Pack × 5 Pcs.)
Dimensions (W × D × H)	140 mm × 155 mm × 55 mm
Dimensions (W × D)	140 mm × 155 mm
Autoclavable	Yes
Material	Polypropylene
Places	100

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de

















More informations about StarStore 100 Storage Box

100-place polypropylene storage box with fixed, 10x10 dividers for tubes up to 13mm, eg. 1.5 / 2.0 ml microcentrifuge tubes or screw cap microtubes. Hinged, lockable lid with positioning lugs for stable stacking. Can be used in the freezer racking made specifically for use with this storage box.



All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	StarStore 100-Place Storage Box, Natural Color: Natural	5 Pcs. (1 Pack × 5 Pcs.)	I2310-5800
	 Create PDF Data Sheet		
	StarStore 100-Place Storage Box, Neon Blue Color: Neon blue	5 Pcs. (1 Pack × 5 Pcs.)	I2310-5841
	 Create PDF Data Sheet		
	StarStore 100-Place Storage Box, Neon Green Color: Neon green	5 Pcs. (1 Pack × 5 Pcs.)	I2310-5842
	 Create PDF Data Sheet		
	StarStore 100-Place Storage Box, Neon Orange Color: Neon orange	5 Pcs. (1 Pack × 5 Pcs.)	I2310-5843
	 Create PDF Data Sheet		
	StarStore 100-Place Storage Box, Neon Pink Color: Neon pink	5 Pcs. (1 Pack × 5 Pcs.)	I2310-5844
	 Create PDF Data Sheet		
	StarStore 100-Place Storage Box, Neon Yellow Color: Neon yellow	5 Pcs. (1 Pack × 5 Pcs.)	I2310-5846
	 Create PDF Data Sheet		
	StarStore 100-Place Storage Box, Neon Mixed* Color: Neon mixed	5 Pcs. (1 Pack × 5 Pcs.)	I2310-5848
	 Create PDF Data Sheet		

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de