



StarStore 50 (hinged)



Product Highlights

- 50-place polypropylene storage box with fixed 10x10 dividers for tubes up to 13mm
- Hinged, lockable lid with positioning lugs for stable stacking
- Moulded reference grid on base.
- Temperature use -90°C to 121°C
- Autoclavable
- Dimensions (WxDxH): 140 x 80 x 55 mm

Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	See variations
Pack Size	5 Pcs. (5 Pcs.)
Dimensions (W × D × H)	140 mm × 80 mm × 55 mm
Dimensions (W × D)	140 mm × 80 mm
Max. tube length	13 mm
Autoclavable	Yes
Material	Polypropylene
Color	Various
Operating temperature	-90°C to 121°C
Places	50
For use with	for tubes up to 13mm

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about StarStore 50 (hinged)















50-place polypropylene storage box with fixed, 5x10 dividers for tubes up to 13mm, eg. 1.5/2.0ml microcentrifuge tubes or screw cap microtubes. Hinged, lockable lid with positioning lugs for stable stacking.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 StarStore 50 Neon (Hinged), Natural Color: Natural  Create PDF Data Sheet	5 Pcs. (1 Pack × 5 Pcs.)	I2350-5800
 StarStore 50 Neon (Hinged), Neon Blue Color: Neon blue  Create PDF Data Sheet	5 Pcs. (1 Pack × 5 Pcs.)	I2350-5841
 StarStore 50 Neon (Hinged), Neon Green Color: Neon green  Create PDF Data Sheet	5 Pcs. (1 Pack × 5 Pcs.)	I2350-5842
 StarStore 50 Neon (Hinged), Neon Orange Color: Neon orange  Create PDF Data Sheet	5 Pcs. (1 Pack × 5 Pcs.)	I2350-5843
 StarStore 50 Neon (Hinged), Neon Pink Color: Neon pink  Create PDF Data Sheet	5 Pcs. (1 Pack × 5 Pcs.)	I2350-5844
 StarStore 50 Neon (Hinged), Neon Yellow Color: Neon yellow  Create PDF Data Sheet	5 Pcs. (1 Pack × 5 Pcs.)	I2350-5846
 StarStore 50 Neon (Hinged), Neon Mixed Color: Neon mixed  Create PDF Data Sheet	5 Pcs. (1 Pack × 5 Pcs.)	I2350-5848

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de