



Centrifuge Tube Storage and Transport Box



Product Highlights

- Sturdy polypropylene storage / transport box for centrifuge tubes
- One holds 36 x 15ml tubes; the other holds 16 x 50ml tubes
- Transparent, lift-off lid
- Base can be removed for use in waterbath
- Temperature use -90°C to 121°C
- Autoclavable
- Dimensions (WxDxH): 145 x 145 x 127 mm

Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	See variations
Dimensions (W × D × H)	145 mm × 145 mm × 127 mm
Dimensions (W × D)	145 mm × 145 mm
Autoclavable	Yes
Pack Size	2 Pcs. (2 Pcs.)
Rack Type	Tube Rack
Material	Polypropylene

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about Centrifuge Tube Storage and Transport Box

Two polypropylene storage / transport boxes for centrifuge tubes. One holds 36 x 15ml tubes; the other holds 16 x 50ml tubes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 Centrifuge Tube Storage / Transport Box for 50ml Tubes, Blue Places: 16  Create PDF Data Sheet	2 Pcs. (1 Box × 2 Pcs.)	E9025-1650
 Centrifuge Tube Storage / Transport Box for 15ml Tubes, Green Places: 36  Create PDF Data Sheet	2 Pcs. (1 Box × 2 Pcs.)	E9025-3615

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de