



StarRack 80



Product Highlights

- 80-place polypropylene fraction collection rack with a 5x16 array for tubes up to 11mm, eg. 1.5/2.0ml microcentrifuge tubes
- Alpha numeric reference grid
- The rack does not have a lid, but an optional, natural polypropylene storage box with a hinged lid into which the rack fits can be used for extra security during freezer storage
- Temperature use: -90 °C to 121 °C
- Autoclavable
- Dimensions (WxDxH): 226 x 67 x 28 mm

Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	See variations
Dimensions (W × D × H)	226 mm × 67 mm × 28 mm
Dimensions (W × D)	226 mm × 67 mm
Autoclavable	Yes
Rack Type	Multi-Tube Rack
Material	Polypropylene
Color	Various
Places	80

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about StarRack 80

One of the most popular racks in the lab. Polypropylene provides high chemical resistance. Temperature use: -90 °C to 121 °C. Autoclavable.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	StarRack 80 polypropylene fraction collection rack, Natural Color: Natural	5 Pcs. (1 Pack × 5 Pcs.)	E2380-5000
 Create PDF Data Sheet			
	StarRack 80 polypropylene fraction collection rack, Mixed Color: Mixed	5 Pcs. (1 Pack × 5 Pcs.)	E2380-5008
 Create PDF Data Sheet			
	StarRack 80 polypropylene fraction collection rack, Neon Mixed Color: Neon mixed	5 Pcs. (1 Pack × 5 Pcs.)	E2380-5048
 Create PDF Data Sheet			

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de