




Wire Storage Racks



Product Highlights

- Chemically-resistant epoxy-coated wire racks
- Suitable for use in the hood
- Temperature use: -80 °C to 121 °C
- Autoclavable

Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	See variations
Autoclavable	Yes
Pack Size	1 Piece (1 Piece)
Color	Blue
Operating temperature	-80°C to 121°C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about Wire Storage Racks

Chemically-resistant epoxy-coated wire racks, suitable for use in the hood. Temperature use: -80 °C to 121 °C.
Autoclavable.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 Racks for 15ml centrifuge tubes (15-place) Tube Style: 15 ml Centrifuge Places: 15  Create PDF Data Sheet	1 Piece (1 Rack × 1 Piece)	E2396-1515
 Racks for 50 ml centrifuge tubes (8-place) Tube Style: 50 ml Centrifuge Places: 8  Create PDF Data Sheet	1 Piece (1 Rack × 1 Piece)	E2396-1550

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de