



Wire Storage Racks



Product Highlights

- Chemically-resistant epoxy-coated wire racks
- Suitable for use in the hood
- Temperature use: -80 °C to 121 °C
- Autoclavable

Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH Neuer Höltigbaum 38

22143 Hamburg
Email: info@starlab.de







General Data

Art. No.	See variations
Pack Size	1 Piece (1 Piece)
Autoclavable	Yes
Color	Blue
Operating temperature	-80°C to 121°C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH Neuer Höltigbaum 38

22143 Hamburg Email: info@starlab.de







More informations about Wire Storage Racks

Chemically-resistant epoxy-coated wire racks, suitable for use in the hood. Temperature use: -80 $^{\circ}$ C to 121 $^{\circ}$ C. Autoclavable.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de







All Variations

PRODUCT NAME	:	PACKAGING SIZE	ART. NO.
	Racks for 15ml centrifuge tubes (15-place) Tube Style : 15 ml Centrifuge Places : 15	1 Piece (1 Rack × 1 Piece)	E2396-1515
	Create PDF Data Sheet		
	Racks for 50 ml centrifuge tubes (8-place) Tube Style: 50 ml Centrifuge Places: 8	1 Piece (1 Rack × 1 Piece)	E2396-1550
	Create PDF Data Sheet		

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH

Neuer Höltigbaum 38 22143 Hamburg Email: info@starlab.de

